

HOSPITAL-LEVEL THERAPY NOW AVAILABLE AT HOME

Your doctor has probably talked with you about the importance of airway clearance therapy for maintaining your health and minimizing hospitalizations.

That's the idea behind The Volara™ System. The Volara System uses Oscillation and Lung Expansion (OLE) technology, which is highly effective for airway clearance, but was never before available outside of hospital settings.

Now, with The Volara System, you can benefit from OLE therapy at home.



ABOUT HILLROM

Hillrom is a global medical technology leader whose 10,000 employees have a single purpose: enhancing outcomes for patients and their caregivers by advancing connected care. Around the world, our innovations touch over 7 million patients each day. They help enable earlier diagnosis and treatment, optimize surgical efficiency and accelerate patient recovery while simplifying clinical communication and shifting care closer to home. We make these outcomes possible through connected smart beds, patient lifts, patient assessment and monitoring technologies, caregiver collaboration tools, respiratory care devices, advanced operating room equipment and more, delivering actionable, real-time insights at the point of care. **Learn more at respiratorycare.hill-rom.com.**

For more information about products or services, please contact Hillrom customer service at 1-800-426-4224.




¹ Huynh TT, Liesching TN, Cereda M, Lei Y, Frazer MJ, Nahouraii MR, Diette GB, Efficacy of Oscillation and Lung Expansion in Reducing Postoperative Pulmonary Complication, Journal of the American College of Surgeons (2019)

Volara™ is a trademark of Hill-Rom Services PTE Ltd. Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.


© 2020 Hill-Rom Services PTE Ltd. ALL RIGHTS RESERVED
APR61302 rev 2 07-MAY-2020 ENG - US



The Volara™ System OLE Therapy



BREATHE EASY.



3 PROVEN THERAPIES¹ 1 VERSATILE DEVICE

The Volara™ System combines three key therapies in a single device. These therapies work together to expand your lungs, mobilizing mucus from your airways.

- **CPEP**
Continuous Positive Expiratory Pressure expands the lungs and holds them open.
- **CHFO**
Continuous High Frequency Oscillation loosens and moves mucus from the small airways to the large airways.
- **NEBULIZER**
Delivers medication that helps loosen mucus.

Your doctor will determine what treatment cycle is best for you. Then The Volara System can be programmed, so it's easy to do your therapy exactly as prescribed.



PERSONALIZED PRESETS

The Volara System will be preset by your medical team, for therapy precisely tailored for you.

PORTABLE AND LIGHTWEIGHT

The Volara System is light and battery operated, so it can be used throughout the home and beyond.

INTUITIVE INTERFACE

The Volara System interface is designed for ease of use.

WORKS FOR YOU

The Volara System can be used with a mouthpiece, face mask, tracheostomy or in-line on a ventilator.



HILLROM IS HERE FOR YOU.

Your health is our top priority. So we do everything possible to make sure you're getting the therapy you need, feel comfortable using it, and can get answers to questions whenever they come up.

- Working with insurance providers
- Training in your home
- 24/7 customer service 1-800-426-4224