

# DEPENDABLE BARIATRIC CARE



**Hillrom™**

Tri-Flex II™ Bariatric Bed



# RELIABLE AND DEPENDABLE CHOICE FOR THE CARE OF BARIATRIC PATIENTS

Hospitals around the country face many challenges, including cutbacks, staffing shortages, and limited resources. Because of the unique medical demands, bariatric patients often require extra care, which can place even more strain on resources and staff.



## DEPENDABLE CARE

- Expandable frame (37"- 48") for patient comfort.
- Integrated scale.
- Battery backup eliminates need to find outlet when repositioning or transferring patients.
- Easy to use caregiver controls.

## RELIABILITY

- Low air loss and alternating pressure surface (Synergy® Air Elite) helps prevent and treat pressure injuries.
- Low deck height of 14" and drop down rails help make the transition to and from bed easier.

## PATIENT COMFORT

- Patient-helper trapeze helps patients assist caregiver during repositioning or to help make themselves more comfortable.
- 1,000 lb patient weight capacity.



## TECHNICAL SPECIFICATIONS

|  |   |
|--|---|
| <b>Bed</b>                                       |   |
| Low Position - Sleep Deck to Floor               | 14" (36 cm)   |
| High Position - Sleep Deck to Floor              | 28" (71 cm)   |
| Maximum Overall Width                            | 40.75" (109 cm)   |
| Maximum Overall Length                           | 94" (939 cm)  |
| Patient Weight Capacity                          | 1,000 lbs (453.6 kgs)   |
| <b>Bed Angles</b>                                |   |
| Head Section, Knee Section                       | 0-60°, 0-24°  |
| Max Trend  | 12°   |
| <b>Sleep Surfaces</b>                            |   |
| Sleep Deck Width                                 | 37"-48" (Rental) (93.98 cm-121.92 cm)<br>37"-54" (Capital) (93.98 cm-137.16 cm)   |
| Sleep Deck Length                                | 86" (137.16 cm)   |
| Foam Surface                                     | 86"L x 37"-48"W x 9.4"H (218.44 cm x 93.98-121.92 cm x 23.876 cm)<br>[includes two 5.5" (13.97 cm interchangeable side bolsters)] |
| Low Air Loss Surface<br>(Synergy® Elite Surface) | 86"L x 39" or 48"W x 10"H (218.44 cm x 99.06 cm or 121.92 cm x 25.4 cm)   |
| <b>Caster &amp; Braking System</b>               |   |
| Caster Diameter                                  | 6" (15.26 cm)   |
| Four Wheel Braking                               | No  |
| Steer and Lock                                   | Yes   |
| <b>Scale System</b>                              |   |
| Accuracy   | 1% of patient weight  |
| Scale Weight Capacity                            | 1,000 lbs (453.6 kgs)   |
| <b>Safety Listing</b>                            | UL 60601-1  |
| <b>Power Requirements</b>                        | 120v AC , 60H   |

## STANDARD FEATURES

Expandable Width Frame and Foam Surface  
In-Bed Scale (standard on rental)  
Pendant  
Battery Backup

## OPTIONAL FEATURES

Full Frame Trapeze  
Synergy® Air Elite Surface



**Hillrom™**

For more information, please contact your local Hillrom sales representative at 1-800-445-3730 or order online at [rental.hillrom.com](http://rental.hillrom.com).

**hillrom.com**

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2021 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. 186703 rev 3 02-MAR-2021 ENG – US