



Progressa® Patient Bed

Item Description: Electric Specialty Patient Bed

Model Number: P7500 / PRO-105

JSN Number: M7007

Warranty: 1 Year From Delivery Date

Installation: _____

Proposal Number: _____

Facility: _____

City, State: _____

Estimated Delivery: _____

Delivery Cost: _____

Taxes: _____

Product Standard Features

SlideGuard® Frame Technology

Patient Controls Backlighting, Night Light

FlexAfoot™ Retractable Foot Control

Advanced MicroClimate® Technology (*Permanently Installed*)
– Standard with Pulmonary and Therapy Surfaces

30°/45° Head-of-Bed Alarm

Bed Controls, In-Bed Scale, Line Managers

3-Mode Bed Exit Alarm (PPM)

Caregiver Pendant

Standard Head and Foot Panel

Foot Controls Module

Oxygen Tank Holder Module

Steering Caster, Accessory Outlet

HandsFree® Emergency CPR, Trendelenburg

NaviCare® Ready – Local Alarm On/Off

Point-of-Care® Siderails Control

Both Siderails with Point-of-Care® Touch Screen

Smart Bed Ready (SideCom® System Required)

6" Urethane Central Locking Caster System

Obstacle Detect System at Both Sides of the Bed

Maximum Patient Weight: 500 lbs

Drainage Bag Holders

Overall Weight: 800 lbs, Safe Working Load: 650 lbs

Four IV Pole Sockets

Frame/Surface Options

Pulmonary Surface (Powered Air)

Therapy Surface (Powered Air)

Prevention Surface (Non-Powered Air)

Frame Only

Chair Egress Frame Design

Dining Chair Frame Design

StayInPlace™ Patient Migration Management System

Rotation (Permanently Installed)

Percussion/Vibration (Permanently Installed)

Options

NUL-Siderail Nurse Call, Universal TV and Lighting Controls

Siderail Nurse Call

Bed Exit alarm to Nurse Call System

Fifth Wheel Steering

IntelliDrive® Powered Transport Mechanism

Progressa® Transport Shelf

Permanent IV Pole (Choose Right or Left Side)

Wireless Connectivity


Progressa® Patient Bed

Item Description: Electric Specialty Patient Bed

Model Number: P7500 / PRO-105

JSN Number: M7007

Proposal Number: _____



Technical Specifications

Length	
Foot Section Retracted	88" (223.5 cm)
Foot Section Extended	98" (249 cm)
Width	40.5" (103 cm)
Highest Position: Top of Seat Section to Floor	35.7" (90.7 cm)
Lowest Position: Top of Seat Section to Floor	16.5" (41.9 cm)
Standard Caster	6" (15.24 cm)

Mattress Details

Prevention Surface	
Mattress Width	35" (88.9 cm)
Mattress Length	84" (213 cm)
Mattress Height	7.1" (18.1 cm)
Mattress Weight	31 lbs (14.1 kgs)

Safe Working Load	650 lbs (295 kgs)
Maximum Patient Weight	500 lbs (227 kgs)
C-Arm Clearance	7-15" (17.8-38.1 cm)
Electrical Information	100/240 V 6/12A 50/60 Hz
Fluid Ingress Protection	IPX4
Radiolucent Sleep Deck	Std 17.7" x 23" (45 cm x 58.5 cm)
Place of Manufacture	USA

Therapy Surface

Mattress Width	35.5" (90.2 cm)
Mattress Length	84" (213 cm)
Mattress Height	8.34" (21.2 cm)
Mattress Weight	45 lbs (20.5 kgs)

Articulation Specifications

Head	0°-67°
Thigh	0°-30°
Foot	0°-70°
Trendelenburg	13°
Tilt Table Function (Reverse Trendelenburg)	-18°

Pulmonary Surface

Mattress Width	35.5" (90.2 cm)
Mattress Length	84" (213 cm)
Mattress Height	8.34" (21.2 cm)
Mattress Weight	49 lbs (22.3 kgs)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2018 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
202063 rev 2 21-MAY-2018 ENG - US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com

Enhancing outcomes for patients and their caregivers:

Hill-Rom