



TS 7000DV SURGICAL TABLE

Item Description: Mobile OR Table

Model Number: 7000dV

JSN Number: _____

Installation: _____

Proposal Number: _____

Facility: _____

City, State: _____

Estimated Delivery: _____

Delivery Cost: _____

Taxes: _____

TABLE PACKAGE

TS 7000dV Table Package includes Table, Double Joint Head Section with pad, One-Part Leg Section with pad, 2 Armboards with pads, Lightweight Transfer Board with pad, Component Cart, 1 Patient Restraint Strap, 1 Rotatable Footboard with pad and 2 clamps, dV Remote Control, Charging Cable and dV Connection Cable, Software for Integrated Table Motion, and 2 Year Service Coverage including Software Upgrade(s) providing wireless capability for use with the *da Vinci Xi*® Surgical System.

OPTIONAL COMPONENTS & ACCESSORIES

- Two-Part Leg Section (Split Leg) with pad
- Split Leg Footboards with pads
- Pelvic Extension with pad
- Trendelenburg Restraint
- Adjustable Footboard with pad
- Lithotomy Leg Stirrups

PRODUCT STANDARD FEATURES

Integrated Table Motion enables the patient to be dynamically positioned while the surgeon operates

Specialized remote guides the OR team through setup and intraoperative use

Wireless connectivity between the surgical system and table

1,000 lb patient weight capacity in standard configuration, including slide

Radiolucent tabletop and slide range for today body intraoperative imaging

Intuitive touchscreen Remote Control and Column Control Panel

Extreme height range for versatile positioning

Self-aligning Hook Coupling system for enhanced safety in attaching, replacing, and removing components

Synchronized leveling system safely and anatomically returns all articulated components to level position

5th Wheel Assist for ergonomic, directional steering of the table in transport and rotation

Overload protection with audible and visual positioning warnings on remote control for enhanced safety

Collision protection help prevent the tabletop from colliding with the floor or base of the table, even in extreme positions

Integrated X-ray Cassette Channel for full body fixed imaging

Double-layer Visco-elastic pads help to reduce pressure in vulnerable areas

Stainless steel, Y-shaped base accommodates foot access and optimal patient access



For more information, please contact your local distributor or Hillrom sales representative.

**TS 7000DV
SURGICAL TABLE**

Item Description:
Mobile OR Table

Model Number:
7000dV

JSN Number:



TECHNICAL SPECIFICATIONS

Maximum Weight Capacity	1,000 lb (450 kg)
Maximum Weight Capacity Fully Articulated	1,000 lb (450 kg)
Height Adjustment	23" - 45" (± 0.4") (59 cm - 114 cm (± 1 cm))
Lift	22" (± 0.4") (55 cm (± 1 cm))
Trendelenburg / Reverse Trendelenburg	± 45°
Tilt	± 30°
Longitudinal Slide	18" (± 4") 11" to the foot; 7" to the head (46 cm (± 1 cm))
Adjustment Range of Leg Section Joints	+ 90° / - 105°
Adjustment of Lower Back Section Joints	+ 90° / - 45°
Warranty	Two-year warranty included with value added SmartCare Protection Plus
Electrical Information	100-230 V, 50/60 Hz 240 V, 50Hz



hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.