



TS 7500 OR SYSTEM TABLES

Item Description: OR System Table

Model Number: TS7500

JSN Number: _____

Installation: _____

Proposal Number: _____

Facility: _____

City, State: _____

Estimated Delivery: _____

Delivery Cost: _____

Taxes: _____

PRODUCT STANDARD FEATURES

A variety of tabletop options provide the user with the right tool for the right case.

Intuitive user interfaces for easy operation in daily surgical work.

Lightweight components and intuitive coupling points require little effort when setting up and configuring the table.

Using force sensor technology, collision and overload protection help to enhance patient safety.

Innovative solutions for hybrid surgery that communication with CT, MRI and angiography systems.

TABLE COLUMN OPTIONS

- TS 7500 System Table Mobile Column
- TS 7500 System Table Fixed Column (in-floor install)
- TS 7500 System Table Fixed Column (on-floor install)

AVAILABLE TABLETOPS

- Universal Tabletop (U14H)
- Carbon FloatLine Tabletop
- Universal Tabletop (U24H)
- MR Neuro Tabletop
- Universal Tabletop (U26H)
- Specialty Tabletop (ST26H)
- Pediatric Tabletop
- Carbon Spine Tabletop

TECHNICAL SPECIFICATIONS

Versions	TS 7500 Mobile Column	TS 7500 Fixed Column (in-floor)	TS 7500 Fixed Column (on-floor)
Maximum Capacity	Static: 880 lbs (400 kg) Articulated: Dependent on tabletop configuration		
Vertical Range	21.7" (55 cm) total lift 20.1" - 41.7" (51-106 cm)	21.7 (55 cm) total lift 19.3" - 40.9" (49-104 cm)	21.7 (55 cm) total lift 20.5" - 42.1" (52-107 cm)
Trendelenburg / Reverse Trendelenburg	± 45° (± 65° with ST26 Tabletop)		
Tilt	± 30°		
Warranty	One-year		
Electrical Information	100-230 V, 50/60 Hz 240 V, 50Hz		



For more information, please contact your local distributor or Hillrom sales representative.

TS 7500 OR SYSTEM TABLE

Item Description: OR System Table

Model Number: TS7500

JSN Number: _____

TECHNICAL SPECIFICATIONS - TABLETOPS

Tabletop	U14H	U24H	U26H	Pediatric	Carbon FloatLine	MR Neuro	ST26H	Carbon Spine
Max. Capacity - Static	880 lbs (400 kg)			176.37 lbs (80 kg) without extension unit; 297 lbs (135 kg) with extension unit & OR Adapter	496 lbs (225 kg)	396.8 lbs (180 kg)	792 lbs (360 kg)	496 lbs (225 kg)
Max. Capacity - Articulated	880 lbs (400 kg) over center and up to 15 degrees tilt; 550 lbs (250 kgs) without restrictions			176.37 lbs (80 kg)	396 lbs (180 kg) without restrictions		792 lbs (360 kg) over center and up to 15° tilt; 550 lbs (250 kg) without restrictions	396 lbs (180 kg)
Trendelenburg / Reverse Trendelenburg	± 45°						± 65°	± 20°
Tilt	± 30°							± 15°
Vertical Range	Mobile: 24.09" - 45.75" (612 - 1162 mm); Fixed (in-floor): 23.31" - 44.96" (592 - 1142 mm); Fixed (on-floor): 24.49" - 46.14" (622 - 1172 mm)				Mobile: 27.52" - 49.17" (699 - 1249 mm); Fixed (in-floor): 26.73" - 48.39" (679 - 1229 mm); Fixed (on-floor): 27.91" - 49.57" (709 - 1259 mm)		Mobile: 26.22" - 47.87" (666 - 1216 mm); Fixed (in-floor): 27.0" - 48.66" (686 - 1236 mm); Fixed (on-floor): 27.40" - 49.05" (696 - 1249 mm)	
Back Section	+90 to -105	+90 to -55 lower back	+90 to -50 lower back +90 to -105 upper back	N/A		+30	+90 to -50 lower back +90 to -105 upper back	N/A
Foot / Leg Section	+90 to -105			+90 / - 60	N/A		+90 to -105	N/A
Head Section	+90 to -45			N/A			+90 to -45	Dependent on head section
Accessory Siderail	Yes			Removeable siderails			Varied	No siderails
Perineal Cutout	Yes			N/A			Yes	N/A
Radiolucent Tabletop	Yes, dependent on configuration; tabletop slide of 15.75" (400 mm).			Yes; 600 mm to the head and foot.	Yes; 82.7" (210 cm) of 360 imaging w/ head positioning accessory	Yes, 360 tabletop rotation	N/A	Yes, 360 tabletop rotation
C-arm Accessible	Yes							
Cassette Sizes	N/A							
Kidney Elevator	No, body elevator			N/A			No, body elevator	N/A
Sliding Top	15.75" (400 mm)			23.62" (600 mm) to head and foot	Mobile: 25.6" (650 mm); Fixed: 31.4" (80 cm); plus Transverse: 9.8" (25 cm)	18.6" (47.2 cm) (+/- 5" (125 mm))	15.75" (40 cm)	Mobile: 23.6" (600 mm); Fixed: 31.5" (800 mm)



hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2020 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. APR117801 rev 1 04-AUG-2020 ENG - US