



TS 7000 OPERATING TABLE

Item Description: Mobile OR Table

Model Number: 7000 / 7000 MB

JSN Number: _____

Installation: _____

Proposal Number: _____

Facility: _____

City, State: _____

Estimated Delivery: _____

Delivery Cost: _____

Taxes: _____

PRODUCT STANDARD FEATURES

1,000 lb patient weight capacity in standard configuration, including slide

Radiolucent tabletop and slide range for total body intraoperative imaging

Intuitive touchscreen Remote Control and Column Control Panel

Extreme height range for versatile positioning

Self-aligning Hook Coupling system for enhanced safety in attaching, replacing, and removing components

Synchronized leveling system safely and anatomically returns all articulated components to level position

5th Wheel Assist for ergonomic, directional steering of the table in transport and rotation

Overload protection with audible and visual positioning warnings on remote control for enhanced safety

Collision protection help prevent the tabletop from colliding with the floor or base of the table, even in extreme positions

Integrated X-ray Cassette Channel for full body fixed imaging

Double-layer Visco-elastic pads help to reduce pressure in vulnerable areas

Stainless steel, Y-shaped base accommodates foot access and optimal patient access

CONFIGURATIONS

TS 7000 Table

TS 7000 Motorized Base (MB) Table

ACCESSORIES

Neurosurgery

Spine Surgery

Orthopedic Surgery

Urology & Gynecology

Bariatric Surgery

Pediatric Applications

Anesthesia

Radiology



For more information, please contact your local distributor or Hillrom sales representative.

**TS 7000
SURGICAL TABLE**

Item Description:
Mobile OR Table

Model Number:
7000

JSN Number:



TECHNICAL SPECIFICATIONS

	TS 7000 / 7000 MB
Maximum Weight Capacity	1,000 lb (450 kg)
Maximum Weight Capacity Fully Articulated	1,000 lb (450 kg)
Height Adjustment	23" - 45" (± 0.4") (59 cm - 114 cm (± 1 cm))
Lift	22" (± 0.4") (55 cm (± 1 cm))
Trendelenburg / Reverse Trendelenburg	± 45°
Tilt	± 30°
Longitudinal Slide	18" (± 4") 11" to the foot; 7" to the head (46 cm (± 1 cm))
Adjustment Range of Leg Section Joints	+ 90° / - 105°
Adjustment of Lower Back Section Joints	+ 90° / - 45°
Warranty	One-year
Electrical Information	100-230 V, 50/60 Hz 240 V, 50Hz



hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.