



Hillrom™



2021 Product Catalog

HEALTHCARE FURNITURE
Distributed by Hillrom

Recliner

FEATURES AND BENEFITS:

- Infinite position recline
- Independent back and ottoman recline
- Robust flip-up arms with locking in all positions up/down
- Side retractable accessory hook
- Caregiver ergonomic push bar on back
- Front and rear lock casters with toe lock/release casters in front
- Compatible with Hillrom™ overbed tables and patient lifts
- Structure composed of injection molded plastic and cast aluminum with no exposed wood
- Molded foam cushions
- Solid surface arm caps available in multiple colors
- Available in non-porous fabrics with minimal sewing seams
- Optional central lock and steer

Dimensions	W: 33.75" D: 35" H: 44.25" SW: 23" SD: 20" SH: 19" Arm Height: 26.25" Back Seat Height: 28" Floor Clearance: 4.5"
Weight Capacity	500 lbs
Caster Size	4"
Range of full recline	74°
Part Number	Non-Powered: P9180 Powered: P9181
Overall Weight	Non-Powered: 114 lbs Powered: 125 lbs



Powered Recliner

- Independent power recliner and ottoman operation, with home return feature
- One integrated USB port
- Battery back-up



Infinite position recline

Convenient, knee-height arm lock push button allows for hands-free engagement for caregivers

Bariatric Recliner

FEATURES AND BENEFITS:

- Infinite position recline up to 45°
- Independent back and ottoman recline
- Side retractable accessory hook
- Caregiver ergonomic push bar on back with release button
- Front and rear lock casters with toe lock/release
- Compatible with Hillrom™ overbed tables and patient lifts
- Structure composed of steel and hardwood with no exposed wood
- Molded foam cushions
- Solid surface arm caps available in multiple colors
- Available in non-porous fabrics with minimal sewing seams
- Optional removable arm (left side only)



Dimensions	W: 36" D: 36.5" H: 45" SW: 29.75" SD: 20" SH: 20.5" Arm Height: 26.25" Back Seat Height: 30"
Weight Capacity	880 lbs
Caster Size	4"
Range of full recline	45°
Part Number	Non-Powered: P9190 Powered: P9191
Overall Weight	Non-Powered: 175 lbs Powered: 185 lbs

Powered Recliner

- Independent power recliner/ottoman
- One integrated USB port
- Battery back-up





Compatible with
Hillrom™ Mobile Patient Lifts



Side retractable
accessory hook

Bariatric Recliner fabric shown is EZ6 Maya

20" Highback and Lowback

FEATURES AND BENEFITS:

- Ergonomic lumbar support back design for comfort
- Non-marring glides
- Wall saver back design
- Structure composed of steel and hardwood with no exposed wood
- Solid surface arm caps and powder coated steel frame available in multiple colors
- Available in non-porous fabrics with minimal sewing seams
- Optional flex back

Dimensions	Lowback: W: 23.5" D: 24.25" H: 33.75" SW: 19.5" SD: 18" SH: 18.5" Arm Height: 25.125"
	Highback: W: 23.5" D: 24.25" H: 40.75" SW: 19.5" SD: 18" SH: 18.5" Arm Height: 25.125"
Weight Capacity	300 lbs
Part Number	P9120
Overall Weight	Lowback: 28 lbs Highback: 30 lbs




Folding Chair with Wall Mount

FEATURES AND BENEFITS:

- Hinged powder coated steel frame available in multiple colors
- Wall mounted hanging bracket matching steel frame color
- Non-marring glides
- Wall saver back design
- Available in seamless, non-porous fabrics
- Mounting hardware for bracket included



Dimensions	Chair: W: 20" D: 19" H: 32.6875" SW: 17.75" SD: 17.75" SH: 17.75" Bracket: W:10" D: 2.15" H: 2.06"
Weight Capacity	300 lbs
Part Number	P9111
Overall Weight	Chair: 22 lbs Bracket: 1 lb



2000+ FABRIC OPTIONS

See page 18

30" and 40" Lowback

FEATURES AND BENEFITS:

- Ergonomic lumbar support back design for comfort
- Non-marring glides
- Wall saver back design
- Structure composed of steel and hardwood with no exposed wood
- Solid surface arm caps and powder coated steel frame available in multiple colors
- Available in non-porous fabrics with minimal sewing seams

30" Lowback

Dimensions	W: 33.5" D: 24.25" H: 33.75" SW: 29.5" SD: 18" SH: 18.5" Arm Height: 25.125"
Weight Capacity	400 lbs
Part Number	P9130
Overall Weight	37 lbs

40" Lowback

Dimensions	W: 43.5" D: 24.25" H: 33.75" SW: 39.5" SD: 18" SH: 18.5" Arm Height: 25.125"
Weight Capacity	500 lbs
Part Number	P9140
Overall Weight	47 lbs





Ergonomic lumbar support
back design for comfort

Structure composed of
steel and hardwood with
no exposed wood

Sleeper Chair

FEATURES AND BENEFITS:

- 2 integrated USB ports
- Easy pull-out mechanism for chaise or sleep surface
- Interior cubby for personal storage while sleeping (24"W x 6"H x 4.625"D)
- Open bottom frame with 4 casters (2 front locking)
- Solid surface arm caps available in multiple colors
- Optional swivel, solid surface tray table (right side sitting)
- Optional LED reading light (right side sitting)

35" Sleeper Chair

Dimensions	W: 35" D: 37" H: 36" SW: 27" SD: 20" SH: 18" Arm Height: 24.5" Chaise: D: 64.125" Sleep Surface: D: 72"
Weight Capacity	350 lbs
Part Number	P9135
Overall Weight	170 lbs




Integrated USB ports



45" Sleeper Chair


Dimensions	W: 45" D: 37" H: 36" SW: 37" SD: 20" SH: 18" Arm Height: 24.5" Chaise: D: 64.125" Sleep Surface: D: 72"
Weight Capacity	350 lbs
Part Number	P9145
Overall Weight	220 lbs



2000+ FABRIC OPTIONS
See page 18



Optional LED
reading light



Interior cubby for
personal storage
while sleeping



Sleeper Sofa

FEATURES AND BENEFITS:

- Available in 76" or 86" width
- 2 integrated USB ports
- Easy fold down sleep pad
- Interior cubby for personal storage while sleeping (24"W x 6"H x 4.625"D)
- Moisture barrier on sleep surface
- Open bottom frame with 5 aluminum legs and non-marring glides
- Solid surface arm caps available in multiple colors
- Optional 5 casters with 2 front locking
- Optional 2 drawers for storage of pillows & blankets (24" W x 6" H x 22" D)
- Optional swivel, solid surface tray table (right side sitting)
- Optional LED reading light (right side sitting)

76" Sleeper Sofa

Dimensions	W: 76" D: 36" H: 33" SW: 67.5" SD: 21.5" SH: 17" Arm Height: 23.75" Sleep Surface: W: 67.5" D: 30"
Weight Capacity	750 lbs
Part Number	P9176
Overall Weight	175 lbs

86" Sleeper Sofa

Dimensions	W: 86" D: 36" H: 33" SW: 77.5" SD: 21.5" SH: 17" Arm Height: 23.75" Sleep Surface: W: 77.5" D: 30"
Weight Capacity	750 lbs
Part Number	P9186
Overall Weight	182 lbs



Integrated USB ports



Interior cubby for personal storage



Dual Top Overbed Table

FEATURES AND BENEFITS:

- Two seamless thermofoil table surfaces:
 - Main table with spill guard edge
 - Flip table with 2 cup holders indentations
- Variable height adjustment from 28.5" to 42.5"
- Integrated rail on both sides of column for optional quick mount accessories
- U-base design for fit and function with beds & furniture
- Low base compatible with beds & furniture
- Graphite grey steel base and steel column with no plastic cover for easy cleaning
- One hand activation to lift flip table
- Available in multiple thermofoil colors

Dimensions	Height: (hi/hi) 42.5" (low/low) 28.5"
	Width: Table: 17" Base: 19"
	Length: Main Table: 32" Flip Table:13.5" Base: 32.5"
	Base Height: 1.8" Table Thickness: 0.75"
Weight Capacity	100 lbs evenly distributed
Part Number	P009498
Overall Weight	55 lbs





Single Top Overbed Table

FEATURES AND BENEFITS:


- Seamless thermofoil table surface with spill guard edge
- Variable height adjustment from 28.5" to 42.5"
- Integrated rail on both sides of column for optional quick mount accessories
- U-base design for fit and function with beds & furniture
- Low base compatible with beds & furniture
- Graphite grey steel base and steel column with no plastic cover for easy cleaning
- Available in multiple thermofoil colors

Dimensions	Height: (hi/hi) 42.5" (low/low) 28.5" Width: Table: 17" Base: 19" Length: Table: 32" Base: 32.5" Base Height: 1.8" Table Thickness: 0.75"
Weight Capacity	100 lbs evenly distributed
Part Number	P009497
Overall Weight	50 lbs

Overbed Table Accessories

Locking Casters

- 2 locking casters at column-side of base




Urinal Hook

- Hook for safe and sanitary urinal storage
- Attaches to integrated column rail
- Note: Disposable urinal not included



Vanity Mirror & Drawer

- Located under table
- Dual-sided for access on either side of table
- Flip up mirror with drawer area for storage
- Drawer Size: 8.8"L x 5.2"W x 2.2"H
- Mirror Size: 9.25"L x 5.5"W



Basket

- Removable storage that allows belongings to be visible
- Attaches to integrated column rail
- 10"L x 6"W x 9"H





Metropolitan®
P528 / BSCM

Dimensions	Metropolitan®: W-22.875", D-16.875", H-31.5" Aero®: W-20", D-17.5", H-30"
Weight Capacity	25 lbs
Part Number	P528-Metropolitan P526-Aero
Overall Weight	92 lbs-Metropolitan 68 lbs-Aero

Art of Care® Bedside Cabinets

FEATURES AND BENEFITS:

- Hardwood and hardwood veneer drawer fronts and side panels for style
- Ventilated back
- Steel frame and base
- Easy roll casters
- Three drawer with easy grip drawer handles
- Ultra tough top, stain and scratch resistant
- Rounded drawer corners for ease of cleaning
- Available in multiple thermofoils and laminates



Aero®
P526 / BSCR

Fabric & Finish Studio

DESIGN YOUR OWN:

This feature gives you the ability to view fabric and finish offerings and "virtually" apply them to the product for easy viewing. You can also order fabrics and finishes at no charge.




FABRIC OPTIONS:

With nearly 2000 options, Hillrom's fabric options are as versatile as the needs of today's healthcare environments. Choose from a wide selection of coated and performance fabrics, ranging from faux leathers to biophilic designs.

HILLROM OFFERS MULTIPLE SOLID COLOR OPTIONS, THERMOFOIL, WOOD, AND POWDER COAT FINISH OPTIONS TO MATCH YOUR HOSPITAL'S DESIGN AND AESTHETICS.





Visit hillrom.com to "Design Your Own" furniture.

2000+
Fabric Options





Hillrom™

ABOUT HILLROM

Hillrom is a global medical technology leader whose 10,000 employees have a single purpose: enhancing outcomes for patients and their caregivers by advancing connected care. Around the world, our innovations touch over 7 million patients each day. They help enable earlier diagnosis and treatment, optimize surgical efficiency and accelerate patient recovery while simplifying clinical communication and shifting care closer to home. We make these outcomes possible through connected smart beds, patient lifts, patient assessment and monitoring technologies, caregiver collaboration tools, respiratory care devices, advanced operating room equipment and more, delivering actionable, real-time insights at the point of care. Learn more at hillrom.com.

For more information, please contact your local distributor or Hillrom sales representative at 1-800-535-6663.

hillrom.com

130 E. Randolph St. Suite 1000, Chicago, IL 60601

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2021 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. APR176801 rev 4 10-DEC-2021 ENG – US