




Hillrom™



Viking® Mobile Lifts
M, L, and XL



Viking M

Landau

Viking M
Lifting Device
Model No. VIK-M-100
Weight Capacity: 300 lbs
Made in USA
© 2010 Viking Medical
Landau

A person with long red hair, wearing blue scrubs, is shown from the back, standing next to a light-colored door. The door has a silver handle and a jagged, sawtooth-shaped cutout on its right side. The background is white.

SAFE PATIENT HANDLING & MOBILITY CALLS FOR A PATIENT- CENTERED APPROACH

Every day, you tailor care and support for each of your unique patients. What if your equipment could do the same?

With the Viking® lifts from Hillrom, you can choose just the right lift for each scenario, and rest assured you'll receive a consistent, high-quality experience.

Available in a variety of sizes and lifting capacities, there is a Viking lift for virtually every need.

UNIQUE SUPPORT FOR EVERY PATIENT

Viking® M Mobile Lift

Compact and designed for your most common lifting needs.

Safe working load: 205 kg / 450 lbs



Viking® L Mobile Lift

A cost-effective solution for larger patient populations.

Safe working load: 250 kg / 550 lbs



EASY TO USE

- Intuitive user display includes lift count, service required, battery usage and overload indicator
- Optional Quick-Release Hook System™ makes it easy to change slingbars or add a scale
- Ergonomic hand control is easy to attach and access from the maneuvering handles
- Two lift speeds support a variety of controlled lifting needs (Viking L and XL only)

LIGHTWEIGHT AND MOBILE

- Lightweight design helps you bring support directly to your patient — when and where they need it
- Optional armrests support patients while they walk
- Large casters roll easily and are easy to clean

DESIGNED FOR SAFETY

- Safe working loads of up to 300 kg / 660 lbs help you lift even your heaviest patients
- Outer Tube™ actuator protection gives added safety
- Accessories like slings, slingbars and scales help you adapt the lift for different needs — and can be used across all three models

Viking® XL Mobile Lift

Powerful enough to support even your largest patients.

Safe working load: 300 kg / 660 lbs



The Support You Need, When and Where You Need It

As your patients move along the mobility continuum and throughout different care settings, their lifting needs change as well. We designed the Viking line of mobile lifts to accommodate virtually all of them — giving you the support you need, every time.

LIFTING TO AND
FROM BED OR CHAIR



TREATMENT AND
REPOSITIONING



LIFTING FROM
THE FLOOR



AMBULATION AND
GAIT TRAINING SUPPORT



VIKING® MOBILE LIFTS TECHNICAL SPECIFICATIONS

Product Number	Viking M Lift 2040045	Viking L Lift 2040044	Viking XL Lift 2040043
Maximum Load	205 kg / 450 lbs	250 kg / 550 lbs	300 kg / 660 lbs
Weight	30 kg / 66 lbs	36.7 kg / 81 lbs	39.9 kg / 88 lbs
Height	1340 - 2020 mm / 53" - 79.5"	1390 - 2030 mm / 54.7" - 80"	1500 - 2110 mm / 59.1" - 83.1"
Base			
Height	106.6 mm / 4.2"	140 mm / 5.5"	140 mm / 5.5"
Length	1230 mm / 48.4"	1360 mm / 53.5"	1420 mm / 55.9"
Width	690 - 1110 mm / 27.2" - 43.3"	695 - 1070 mm / 27.4" - 42.1"	815 - 1130 mm / 32.1" - 44.5"
Expansion	Electric	Electric	Electric
Casters			
Front	75 mm / 3" twin	100 mm / 4" twin 75mm / 3" twin	100 mm / 4" twin 75mm / 3" twin
Rear	75 mm / 3" twin	125 mm / 5" twin 100 mm / 4" twin	125 mm / 5" twin 100 mm / 4" twin
Control Box			
Audible low battery signal	When <25% power remains ^{1,2}	When <25% power remains ^{1,2}	When <25% power remains ^{1,2}
Service Information	Displayed	Displayed	Displayed
Lift counter	Displayed	Displayed	Displayed
Overload indicator	Displayed ²	Displayed ²	Displayed ²
Battery power indicator	Displayed ²	Displayed ²	Displayed ²
Lifting Speed (without load)	36 mm / 1.4" per sec	23 mm / 0.9" per sec and 17mm / 0.7" per sec	23 mm / 0.9" per sec and 17mm / 0.7" per sec
Maximum Lifting Height	1790 mm / 70.5"	1790 mm / 70.5"	1910 mm / 75.2"
Lifting Interval	1270 mm / 50"	1330 mm / 52.4"	1370 mm / 54"
Emergency Lowering	Electric and manual	Electric and manual	Electric and manual
Sound Level	46 dbA	51 dbA	51 dbA
Battery	Lithium ion or Lead-acid gel	Lithium ion or Lead-acid gel	Lithium ion or Lead-acid gel
Charger Options	Built-in and Wall unit	Built-in and Wall unit	Built-in and Wall unit
Hand Control	Wired	Wired	Wired
Protection Class	IPX4	IPX4	IPX4
Material	Aluminum	Aluminum	Aluminum
Slingbar Standard	Universal 450 Slingbar	Universal 450 Slingbar	Universal TwinBar 670
ACCESSORIES			
	LikoScale™ 350, 200 ³ , 400 ³ Slingbar Cover Paddy Quick-Release Hook Armrest Lift Stretcher: OctoStretch™	LikoScale™ 350, 400 ³ Slingbar Cover Paddy Quick-Release Hook Armrest Lift Stretchers: FlexoStretch™ OctoStretch™ LikoStretch 600 IC™	LikoScale™ 350, 400 ³ Slingbar Cover Paddy Quick-Release Hook Armrest Lift Stretchers: FlexoStretch™ OctoStretch™ LikoStretch 600 IC™

¹A buzzer will sound when pressing a push button on the hand control.

²Also displayed on the hand control.

³Available in the United States and Canada only.

IF YOU KNOW LIKO, YOU KNOW HILLROM.

**We've unified our brands to advance connected care.
Visit us at hillrom.com.**



ABOUT HILLROM

Hillrom is a global medical technology leader whose 10,000 employees have a single purpose: enhancing outcomes for patients and their caregivers by advancing connected care. Around the world, our innovations touch over 7 million patients each day. They help enable earlier diagnosis and treatment, optimize surgical efficiency and accelerate patient recovery while simplifying clinical communication and shifting care closer to home. We make these outcomes possible through connected smart beds, patient lifts, patient assessment and monitoring technologies, caregiver collaboration tools, respiratory care devices, advanced operating room equipment and more, delivering actionable, real-time insights at the point of care. Learn more at hillrom.com.

For more information, please contact your Hillrom sales representative at 1-800-445-3730.

hillrom.com

Viking™ mobile lifts are versatile lifts intended mainly for use in healthcare, intensive care and rehabilitation. Viking L and XL lifts are intended for heavier patients.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2021 Liko R&D AB. ALL RIGHTS RESERVED. APR238902-EN-r1 17-JUN-2021 ENG – US