



## Centrella™ Smart+ Bed

**Item Description:** Electric Med-Surg Patient Bed

**Model Number:** CNT155 / CNT255 / CNT355

**JSN Number:** \_\_\_\_\_

**Warranty:** \_\_\_\_\_

**Installation:** \_\_\_\_\_

**Proposal Number:** \_\_\_\_\_

**Facility:** \_\_\_\_\_

**City, State:** \_\_\_\_\_

**Estimated Delivery:** \_\_\_\_\_

**Delivery Cost:** \_\_\_\_\_

**Taxes:** \_\_\_\_\_

### Product Standard Features

Enhanced 3-Mode Bed Exit  
 Motion-Activated Nightlight  
 One-button Patient Stand Assist  
 Advanced SlideGuard® Technology  
 Foley Position Limit  
 Patient HOB Angle Display  
 Trendelenburg/Reverse Trendelenburg  
 Damped Emergency CPR  
 Device Storage Area  
 Integrated Urinal Holders  
 FlexAfoot™ Powered Bed Extension  
 Graphical Caregiver Interface (GCI)® Touchscreen (right side)  
 Bed Exit Silence  
 Brake-Not-Set Alert  
 One Button Chair  
 Scale (kg/lbs)  
 Battery Backup  
 Digital Angle Indicator  
 Four-Corner Brake and Steer  
 Line and Cord Management  
 Four IV Pole Sockets (located at each corner)  
 Integrated Restraint Holders  
 Lockout Controls  
 Head-End Traction  
 Drainage Bag Holder  
 Smart Bed Ready (SideCom® Required)  
 WiFi

### Additional Standard Features for CNT255

Verbal Safety Prompts  
 USB Charging  
 Patient Pendant

### Additional Standard Features for CNT355

Centrella™ max Powered Air Surface  
 IllumiGuide™ Siderail Handgrip  
 SafeView®+ System  
 Second Graphical Caregiver Interface (GCI)® Touchscreen

### CNT155 Product Optional Features

- Patient Pendant
- Second Patient Pendant
- Nurse Call
- Universal TV and Lighting Controls
- 5th Wheel Steering
- Accessory Outlet

### CNT255 Product Optional Features

- IllumiGuide™ Siderail Handgrip
- Comfort Edition 40" Width
- Second Patient Pendant
- Nurse Call
- IntelliDrive® Powered Transport
- Universal TV and Lighting Controls
- SafeView®+ System
- 5th Wheel Steering
- Accessory Outlet

### CNT355 Product Optional Features

- Comfort Edition 40" Width
- Second Patient Pendant
- Nurse Call
- IntelliDrive® Powered Transport
- Universal TV and Lighting Controls
- 5th Wheel Steering
- Accessory Outlet


## Centrella™ Smart+ Bed

**Item Description:** Electric Med-Surg Patient Bed

**Model Number:** CNT155 / CNT255 / CNT355

**JSN Number:** N/A

**Proposal Number:** \_\_\_\_\_



### Technical Specifications

|                                    |  |
|------------------------------------|--|
| Overall width                      | 40" (101.6 cm) or<br>43" (109.22 cm) (Comfort Edition) |
| Overall length with foot extension |  |
| Maximum (foot section extended)    | 98.7" (251 cm)   |
| Minimum (foot section retracted)   | 86.7" (220 cm)   |
| Bed Height maximum                 | 32.5" (83 cm)  |
| Bed Height minimum                 | 14.5" (37 cm)  |
| Head section inclination           | 0° to 65°  |
| Thigh section inclination          | 0° to 25°  |
| Foot section inclination           | 0° to 50°  |

|   |                                     |
|---|-------------------------------------|
| Trendelenburg position (max)                                      | 15°                                 |
| Reverse Trend position (max)                                      | 15°                                 |
| Caster size   type  | 5" (13 cm)   Central Locking System |
| CNT155 core Surface   | 35.5"W x 76"-88"L x 6"H             |
| CNT255 pro Surface  | 35.5" or 38.5"W x 76"-88"L x 7"H    |
| CNT355 max Surface  | 35.5" or 38.5"W x 76"-88"L x 8"H    |
| Safe Working Load   | 650 lbs (295 kg)                    |
| <small>(includes patient weight, surface and accessories)</small> |                                     |
| Patient weight range  | 70 – 500 lbs (32 – 227 kg)          |
| Electrical Information  | 100/240 V   10A   50/60 Hz          |

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2018 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.  
207042 rev 2 21-MAY-2018 ENG – US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

**USA 800-445-3730**  
**Canada 800-267-2337**

[www.hill-rom.com](http://www.hill-rom.com)

Enhancing outcomes for patients and their caregivers:

**Hill-Rom**