



Welch Allyn
Connex® Spot Monitor

Welch Allyn Connex[®] Spot Monitor

Accurate, automatic vitals sent to your Electronic Medical Records (EMR) from the point of care for increased staff productivity, better patient outcomes, and risk reduction.

The only monitor to turn to for truly wireless data transfer, accurate vital signs measurement, and compatibility with your EMR now and in the future.

Accurate readings across your patient population

- > Help ensure accurate hypertension diagnosis with unique SureBP[®] technology with blood pressure averaging.
 - Research shows that more than half of patients could be misdiagnosed in a hypertensive category using the single reading method.¹
- > Provide SpO₂ accuracy in challenging conditions including motion, thick fingers or dark skin tone, and low perfusion with Masimo SET[®] and Nonin PureSAT[®] pulse oximetry technology.²

Improve data accuracy, efficiency, and mobility

- > Eliminate data entry errors in the patient chart with wireless transfer of vital signs.
- > Increase staff efficiency by eliminating typing into patient records.
- > Freedom from wires and cables provides the flexibility to easily work in limited space.





Configurable design

allows you to customise the monitor to the needs of your clinic

Select from these options to get exactly the monitor you need:

Thermometry

- SureTemp® Plus Thermometry
- Braun PRO 6000 Ear Thermometry

Pulse Oximetry

- Masimo SET SpO₂
- Nonin PureSAT SpO₂
- Nellcor® OxiMax™ SpO₂

Mobile Stand

- Accessory Power Management Stand features an integrated battery to ensure the device gets you through a busy day
- Classic Wire Basket Stand (pictured)

EMR partnerships

ensure connectivity over time

- > Welch Allyn closely partners with more than 90 leading EMR providers—including Epic® and Cerner®—to ensure reliable device connectivity.
- > As your EMR evolves, the Connex Spot Monitor will remain compatible with its latest enhancements and changes.

Protect your investment with service programmes

- > **Partners in CareSM Comprehensive Partnership Programmes**—include complete coverage of repair services, protection of key accessories, and annual factory calibration service to help minimise the cost and complexity of servicing and maintaining your monitors.
- > **Unique Enterprise Deployment Solution**—ensures your monitors arrive configured to your specification right out of the box.

¹Robert Smith, MD. [Blood pressure averaging methodology: Decreasing the rate of misdiagnosing hypertension](#). 2014.

²Nonin Medical, [Patient conditions can be unpredictable. Pulse oximeter accuracy shouldn't be](#). 2014

Configurations Available

CONNEX SPOT MONITOR	7100	7300
Parameters		
SureBP Blood Pressure	●	●
Blood Pressure Averaging	●	●
SureTemp Plus Thermometry	○	○
Braun PRO 6000 Ear Thermometry	○	○
Masimo or Nellcor OxiMax SpO ₂		○
Nonin SpO ₂	○	○
Communications		
Wireless (Bluetooth®)		●
USB, Ethernet	●	●
Accessories		
Classic Stand	○	○
Accessory Power Management Stand	○	○

● Included

○ Optional

ACCESSORIES

7000-APM	Connex Spot Accessory Power Management Stand
7000-MS3	Connex Spot Classic Mobile Stand
7000-GCX	Connex Spot Wall Channel
7000-PS	Power Supply (required if not ordering a stand or wall mount)
BATT22	Connex Spot Monitor Lithium-Ion Battery
BATT99	Accessory Power Management Stand Lithium-Ion Battery
4500-925	Connectivity Accessory Kit USB Cable
6000-50	VSM 6000 USB Configuration Memory Stick
4500-35	Blood Pressure Hose with FlexiPort (10 ft)
7000-33	Neonatal Blood Pressure Hose (10 ft)
02895-100	SureTemp Plus Thermometer Rectal Probe and Well
718584	Braun PRO 6000 Ear Thermometer Tether
3278-010	Nonin Reusable Adult Sensor (2 meters)
2360-010	Nonin Reusable Pediatric Sensor (2 meters)
S1-CSM-5	Connex Spot Monitor Comprehensive Partner Program, 5 years*
S2-CSM-5	Connex Spot Monitor Biomed Partner Program, 5 years*

*One- and two-year service plans also available.

Welch Allyn UK Ltd.
Cublington Road Aston Abbotts
Buckinghamshire, HP22 4ND
(Tel:) 0207 365 6780 (Fax) 0207 365 9694

WWW.WELCHALLYN.CO.UK

