





# Hillrom™


HILLROM  
1069 STATE ROUTE 46 EAST  
BATESVILLE, IN 47006  
PHONE: (812) 934-7777  
FAX: (812) 934-8189  
[www.hillrom.com](http://www.hillrom.com)




ISOMETRIC VIEW



PLAN VIEW



FRONT ELEVATION



SIDE ELEVATION


THE VIKING® SERIES OF LIFTS IS READY TO SOLVE MANY VARYING LIFTING NEEDS IN DIFFERENT CARE SITUATIONS, AND THE ADAPTATION POSSIBILITIES ARE AS COMPREHENSIVE AS HILL-ROM'S LIKO® ACCESSORY RANGE. THE VIKING® LIFT SERIES IS AVAILABLE IN FIVE MODELS. FROM THE SMALLEST VIKING® XS LIFT TO THE SUPER-STRONG VIKING® XL LIFT.

THE VIKING® XL LIFT HAS BEEN DEVELOPED FOR VERY HEAVY PATIENTS, BUT IS ALSO SUITABLE FOR AVERAGE WEIGHTED PEOPLE. DESPITE ITS SIZE, THIS MOBILE LIFE WEIGHS ONLY 90 LBS AND CAN LIFT UP TO 660 LBS, MAKING THE VIKING® XL LIFT IDEAL FOR MOST HEAVY-LIFTING SITUATIONS.

© [[2018]] HILLROM SERVICES, INC. ALL RIGHTS RESERVED.

NOTES:

1. INSTALLATION TO BE COMPLETED IN ACCORDANCE WITH MANUFACTURER'S SPECIFICATIONS.
2. DO NOT SCALE DRAWING.
3. THIS DRAWING IS INTENDED FOR USE BY ARCHITECTS, ENGINEERS, CONTRACTORS, CONSULTANTS AND DESIGN PROFESSIONALS FOR PLANNING PURPOSES ONLY. THIS DRAWING MAY NOT BE USED FOR CONSTRUCTION.
4. ALL INFORMATION CONTAINED HEREIN WAS CURRENT AT THE TIME OF DEVELOPMENT BUT MUST BE REVIEWED AND APPROVED BY THE PRODUCT MANUFACTURER TO BE CONSIDERED ACCURATE.
5. CONTRACTOR'S NOTE: FOR PRODUCT AND COMPANY INFORMATION VISIT [www.CADdetails.com/info](http://www.CADdetails.com/info) AND ENTER REFERENCE NUMBER 4794-055



## PATIENT LIFTS (PREVIOUSLY LIKO)

VIKING® XL MOBILE LIFT