



Green Series™ 300 General Exam Light

Item Description: General exam light ideal for patient exam rooms

Model Number: 444XX: Green Series 300 General Exam Light

Installation: _____

Proposal Number: _____

Facility: _____

City, State: _____

Estimated Delivery: _____

Delivery Cost: _____

Taxes: _____

Product Standard Features

A single LED provides whiter, brighter light than halogen

Touchless on/off to help reduce cross-contamination risk

Compact design provides intense light output and cool operation

Less energy consumption and LED life of 50,000 hours helps maximize investment

The Green Series 300 General Exam Light is ideal for: primary care facilities, emergency rooms, urgent care centers, outpatient facilities and dermatology.

Common Configurations

- 44400 Green Series 300 General Exam Light with Mobile Stand
- 44410 Green Series 300 General Exam Light with Table/Wall Mount

Accessories

- 44215 12" Wall Mount Extension
- 48955 Wall/Table Mount
- 48960 Heavy Duty/Tall Mobile Stand (36"/91cm)
- 48950 Mobile Stand (24"/61cm)
- PWCD-B GS 300 Power Cord (8'/2.4m)

#44215
 12" Wall Mount Extension
 (Order separately with #44410 for setup)



24" High
 Standard Stand
 (shipped with #44400)



36" High #48960
 Heavy Duty/Tall
 Mobile Stand
 (Order separately with
 #44410 for setup)



Green Series™ 300 General Exam Light

Item Description: General exam light ideal for patient exam rooms

Model Number:
 444XX: Green Series 300 General Exam Light

Product Specifications

Number of LEDs	1 White LED
Color Temperature	5,500° K
Light Intensity @ Typical Working Distance	12" distance: 30,000 lux (2,787 footcandles) 20" distance: 15,000 lux (1,393 footcandles)
LED Life	50,000 hours
Working Distance	12 in (30.5 cm) to 20 in (50 cm)
Spot Size at Working Distance	5 in (13 cm) to 6 in (15 cm)
Electrical	Input: 100-240 V- 50 Hz-60 Hz, 0.35 A max
Power Cord	#18 AWG, 3-wire grounded cord set, 8' (2.5 m) long
Physical Dimensions	Weight: Total Table/Wall Mount: 5.65 lbs (2.56 kg) Total Mobile Stand: 12.75 lbs (5.78 kg) Control Box: Height: 10.83 in (27.5 cm) Width: 3.94 in (10.0 cm) Depth: 3.15 in (8.0 cm) Luminaire: Height: 7.48 in (19 cm) Width: 3.54 in (9 cm) Depth: 3.15 in (8 cm) Mobile Stand: Height: 24.02 in (61 cm) Diameter: 20.87 in (53 cm) Total Wall Mount: 33.27 in (84.5 cm) Total Mobile Stand: 58.66 in (149.0 cm) Table/Wall and Mobile Stand to Luminaire Center: 31.1 in (79 cm)
Mounting Options	Table/Wall Mount, Mobile Stand
Control Box	Touchless On/Off Control

Welch Allyn, Inc.
 4341 State Street Road
 Skaneateles Falls, NY 13153 USA
 (p) 800.535.6663 (f) 315.685.3361