



Green Series™ 600 Minor Procedure Light

Item Description: Minor procedure exam light designed with many healthcare environments in mind

Model Number: 446XX: Green Series 600 Procedure Light

Installation: _____

Proposal Number: _____

Facility: _____

City, State: _____

Estimated Delivery: _____

Delivery Cost: _____

Taxes: _____

Product Standard Features

Three LEDs that provide a whiter, brighter light than halogen

Less energy consumption for reduced cost of ownership

No bulb replacement for cost savings

Wider head, designed specifically for demanding viewing conditions

Uniform spot and superior brightness for minor procedures and exams in any healthcare environment

Minimal heat dissipation

Touchless controls with dimming features

The Green Series 600 Minor Procedure Light is ideal for: primary care facilities, treatment centers, outpatient facilities, hospital-emergency rooms, dermatology and urgent care centers.

Common Configurations

- 44600 Green Series 600 Minor Procedure Light with Mobile Stand
- 44610 Green Series 600 Minor Procedure Light with Table/Wall Mount

Accessories

- 52630 Disposable Sheath (50/box, 5 boxes/case)
- 44215 12" Wall Mount Extension
- 48960 Heavy Duty/Tall Mobile Stand

#44215
12" Wall Mount Extension
(Order separately with #44410 for setup)



24" High
Standard Stand
(shipped with #44600)

36" High #48960
Heavy Duty/Tall
Mobile Stand
(Order separately with
#44610 for setup)

Green Series™ 600 Minor Procedure Light

Item Description: Minor procedure exam light designed with many healthcare environments in mind

Model Number: 446XX: Green Series 600 Procedure Light

Product Specifications

Number of LEDs	3 White LED
Color Temperature	5,500° K
Light Intensity @ Typical Working Distance	High Intensity Setting: 12" distance: 75,000 lux (6,967 footcandles) 20" distance: 40,000 lux (3,716 footcandles) Low Intensity Setting: 12" distance: 55,000 lux (5,109 footcandles) 20" distance: 30,000 lux (2,787 footcandles)
LED Life	50,000 hours
Working Distance	12 in (30.5 cm) to 20 in (50 cm)
Spot Size at Working Distance	6 in (15 cm) to 7 in (18 cm)
Electrical	Input: 100-240 V- 50 Hz-60 Hz, 0.65 A max
Power Cord	#18 AWG, 3-wire grounded cord set, 8' (2.5 m) long
Physical Dimensions	Weight: Total Table/Wall Mount: 6.45 lbs (2.93 kg) Total Mobile Stand: 17.35 lbs (7.87 kg) Control Box: Height: 10.83 in (27.5 cm) Width: 3.94 in (10.0 cm) Depth: 3.15 in (8.0 cm) Luminaire: Height: 13.58 in (34.5 cm) Width: 7.56 in (19.2 cm) Depth: 3.39 in (8.6 cm) Mobile Stand: Height: 24.02 in (61 cm) Diameter: 20.87 in (53 cm) Total Wall Mount: 33.27 in (84.5 cm) Total Mobile Stand: 58.66 in (149.0 cm) Table/Wall and Mobile Stand to Luminaire Center: 31.1 in (79 cm)
Mounting Options	Table/Wall Mount, Mobile Stand
Control Box	Touchless On/Off Control with Dimming

Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361