



Green Series™ 900 Procedure Light

Item Description: Procedure exam light designed for the most rigorous medical settings.

Model Number: 449XX: Green Series 900 Procedure Light

Installation: _____

Proposal Number: _____

Facility: _____

City, State: _____

Estimated Delivery: _____

Delivery Cost: _____

Taxes: _____

Product Standard Features

- _____ Six LEDs provide whiter, brighter light than halogen
- _____ Touchless on/off and intensity controls help reduce cross-contamination risk
- _____ Extraordinary maneuverability and multiple mounting options
- _____ Uniform spot and superior brightness for demanding viewing conditions
- _____ Less energy consumption and LED life of 50,000 hours helps maximize investment


The Green Series 900 Procedure Light is ideal for: ambulatory care centers, emergency rooms, urgent care centers, intensive care units, labor and delivery and outpatient facilities.

Common Configurations

- 44900 Green Series 900 Procedure Light with Mobile Stand
- 44900-C Green Series 900 Procedure Light with Ceiling Mount
- 44900-W Green Series 900 Procedure Light with Wall Mount

Accessories

- 52630 Disposable Sheath (50/box, 5 boxes/case)
- PWCD-B GS 900 Power Cord (8'/2.4m)



Green Series™ 900 Procedure Light

Item Description: Procedure exam light designed for the most rigorous medical settings.

Model Number: 44900: Green Series 900 Procedure Light

Product Specifications

Number of LEDs	6 White LEDs
Color Temperature	5,500° K
Light Intensity @ Typical Working Distance	<p>High Intensity Setting: 24" distance: 100,000 lux (9,290 footcandles) 40" distance: 55,000 lux (5,109 footcandles)</p> <p>Low Intensity Setting: 24" distance: 75,000 lux (6,967 footcandles) 40" distance: 40,000 lux (3,716 footcandles)</p>
LED Life	50,000 hours
Working Distance	24 in (60.9 cm) to 40 in (100 cm)
Spot Size at Working Distance	9 in (23 cm) to 11 in (28 cm)
Electrical	Input: 100-240 V- 50 Hz-60 Hz, 1.00 A max
Power Cord	#18 AWG, 3-wire grounded cord set, 8' (2.5 m) long
Physical Dimensions	<p>Weight: Total Wall Mount: 17.90 lbs (8.12 kg) Total Mobile Stand: 44.85 lbs (20.34 kg) Total Ceiling Mount: 24.50 lbs (11.11 kg)</p> <p>Luminaire: Height: 13.74 in (34.9 cm) Width: 18.03 in (45.8 cm) Depth: 4.02 in (10.2 cm)</p> <p>Wall to Luminaire Center: 68.9 in (175 cm) Wall Mount Max. Up/Down: 37.4 in (95 cm) Mobile Stand: Height: 68.9 in (175 cm) Mobile Stand Base: Length: 21.97 in (55.8 cm) Width: 22.24 in (56.5 cm) Height: 6.5 in (16.5 cm) Mobile Stand Max. Up/Down: 37.4 in (95 cm) Mobile Stand to Luminaire Center: 37.4 in (95 cm) Ceiling to Luminaire Center: 68.9 in (175 cm) Ceiling Mount Max. Up/Down: 37.4 in (95 cm)</p>
Mounting Options	Wall Mount, Mobile Stand and Ceiling Mount
Control Box	Touchless On/Off Control with Dimming

Welch Allyn, Inc.
 4341 State Street Road
 Skaneateles Falls, NY 13153 USA
 (p) 800.535.6663 (f) 315.685.3361