

Baxter
Hillrom
Progressa+
ICU HOSPITAL BED

**CARE TEAMS
AND PATIENTS
DESERVE
A TRUE ICU BED**



STANDARD FEATURES

- SlideGuard Frame Technology
- Patient Controls Backlighting
- FlexAfoot Retractable Foot Control
- Advanced MicroClimate Technology
- Continuous Low Pressure (CLP)
- Radiolucent Deck
- 30°/45° Head of Bed Alarm
- In-Bed Scale
- Line Managers
- Three Mode Bed Exit Alarm
- Caregiver Pendant
- Caregiver Foot Controls
- Night Light
- Integrated Oxygen Tank Holder
- Accessory Outlet
- Hands Free Emergency CPR & Trendelenburg
- Local Alarm On/Off
- Point of Care Siderails Control
- Both Siderails with Point of Care Touch Screen
- Smart Bed Ready (Sidecom Required)
- Obstacle Detect System at Both Sides of Frame
- Drainage Bag Holders
- Four IV Sockets
- Wireless Connectivity

OPTIONAL FEATURES

- Choice of Support Surface (Essential Surface, Advance Surface, or Accelerate Surface)
- Continuous Lateral Rotation Therapy (CLRT)
- Chest Physiotherapy (CPT)
- Nurse Call, Univ TV & Lighting
- Siderail Nurse Call
- Bed Exit Alarm to Nurse Call
- 5th Wheel Steering
- IntelliDrive Powered Transport Steering
- Transport Shelf
- Permanent IV Pole
- FullChair Position with Sit-to-Stand Option
- StayInPlace Patient Migration Management System
- Alternating Low Pressure (ALP)
- Xtend 40" frame
- Experience POD Device
- USB Charging Option

PROGRESSA+ TRUE ICU BED

TECHNICAL SPECIFICATIONS

| | |
|--|---|
| Total Length | |
| Foot Section Extended | 98" (248.9 cm) |
| Foot Section Retracted | 88" (223.5 cm) |
| Maximum Width - Standard | |
| Side Rails Down | 40.5" (103 cm) |
| Side Rails Up | 40.5" (103 cm) |
| Maximum Width - Xtend | |
| Side Rails Down | 43" (109 cm) |
| Side Rails Up | 43" (109 cm) |
| Maximum Side Rail Height | |
| Above Sleep Deck | 17" (43.2 cm) |
| Lowest Bed Height | |
| Deck to Floor | 16.5" (41.9 cm) |
| Minimum Underbed Clearance | |
| | 4.3" (10.9 cm) w/o IntelliDrive Powered Transport |
| | 1.1" (2.8 cm) w/ IntelliDrive |
| Standard Caster Size | 6" (15.2 cm) |
| Articulations | |
| Head Section | 0° - 67° |
| | 77° with Sit-to-Stand |
| Thigh Section | 0° - 30° |
| Foot Section | 0° - 45° |
| | 75° with Sit-to-Stand |
| Trendelenburg Position (maximum) | 13° |
| Reverse Trendelenburg Position (maximum) | 18° |

| | |
|--------------------------------|--|
| Safe Working Load | 650 lb (295 kg) |
| Maximum Patient Weight | 500 lb (227 kg) |
| C-Arm Clearance | 7 - 15" (17.8-38.1 cm) |
| Cleanability/Infection Control | IPX4 |
| Radiolucent Sleep Deck | Standard 17.7" x 24.8" (45 cm x 63 cm) |
| Voltages | 100/110/115/120/127/220/230/240 |
| Mattress | |
| Essential Surface | |
| Mattress Width | 35" (88.9 cm) |
| Mattress Length | 84" (213 cm) |
| Mattress Height | 7.1" (18.1 cm) |
| Mattress Weight | 31 lb (14.1 kg) |
| Advance Surface | |
| Mattress Width - Standard | 35.5" (90.2 cm) |
| Mattress Width - Xtend | 39" (99 cm) |
| Mattress Length | 84" (213 cm) |
| Mattress Height | 9" (22.86 cm) |
| Mattress Weight | 45 lb (20.5 kg) |
| Accelerate Surface | |
| Mattress Width - Standard | 35.5" (90.2 cm) |
| Mattress Width - Xtend | 39" (99 cm) |
| Mattress Length | 84" (213 cm) |
| Mattress Height | 8" (20.32 cm) |
| Mattress Weight | 48 lb (21.8 kg) |

SURFACE & FRAME CONFIGURATION OPTIONS

| PROGRESSA+ ESSENTIAL | PROGRESSA+ ADVANCE | PROGRESSA+ ACCELERATE |
|---|---|---|
| Essential Surface | Advance Surface | Choose one: - Pulmonary Surface (CLRT only) - Pulmonary Surface (CLRT & CPT) |
| Choose one: - Dining Chair Position - FullChair Position with Sit-to-Stand - FullChair Position with Sit-to-Stand and StayInPlace Technology | Choose one: - Dining Chair Position - FullChair Position with Sit-to-Stand - FullChair Position with Sit-to-Stand and StayInPlace Technology | Choose one: - Dining Chair Position - FullChair Position with Sit-to-Stand - FullChair Position with Sit-to-Stand and StayInPlace Technology |

Baxter.com

Baxter, Advanced MicroClimate, Dining Chair, Experience POD, FlexAFoot, FullChair, Hillrom, IntelliDrive, Obstacle Detect, Point of Care, Progressa+, SlideGuard and StayInPlace are trademarks of Baxter International Inc. or its subsidiaries.