

Hill-Rom® Transport Stretcher

Item Description: Transport Stretcher

Model Number: P8005

JSN Number: M4655

Warranty: 1 Year From Delivery Date

Installation: _____

Proposal Number: _____

Facility: _____

City, State: _____

Estimated Delivery: _____

Delivery Cost: _____

Taxes: _____



Product Standard Features

26" Width

Bumper Protected Siderails

Trendelenburg/Reverse Trendelenburg

Six Transport Strap Locations

3" Comfortline® Foam Mattress

Locking Steering Caster

Easy Release Fold-Down Siderails

Brake/Steer – Head and Foot

Dual-Sided hi/lo Hydraulic Control Pedals

3/4" Standard IV Pole Sockets

Product Optional Features

☐ 30" Width

☐ 4" or 5" Comfortline® Foam Mattress

☐ 5" AccuMax® VPC Surface (Elastimax or Hypolex)

☐ Steering Plus™ System

☐ Brake/Steer – All Four Sides

☐ 3-Sided hi/lo Hydraulic Control Pedals

☐ Ergonomic Push Handles

☐ Ergonomic Push Handles with Integrated IV Pole Transport Holder

☐ Integrated O₂ Tank Holder with Utility Tray

☐ Removable IV Pole – 2 Stage

☐ Permanent IV Pole – 2 Stage

☐ Radiolucent Deck

☐ Upright Chest X-Ray Cassette Holder

Hill-Rom® Transport Stretcher

Item Description: Transport Stretcher

Model Number: P8005

JSN Number: M4655

Proposal Number: _____



Technical Specifications

Maximum weight limit	700 lbs (317.5 kg)
Low position (floor to top of sleep deck)	20.50" (52.07 cm)
High position (floor to top of sleep deck)	34.25" (87 cm)
Overall length	83" (210.8 cm)
Overall width (siderails up)	32" (81.3 cm) or 36" (91.4 cm)
Overall width (siderails stored)	30.38" (77.15 cm) or 34.38" (87.31 cm)
Siderail length	47" (119.38 cm)
Siderail height above sleep deck	14" (35.6 cm)
Comfortline® mattress size	26" (66 cm) x 75" (190.5 cm) or 30" (76.2 cm) x 75" (190.5 cm)
AccuMax® VPC mattress size	26" (66 cm) x 76" (193 cm) or 30" (66 cm) x 76" (193 cm)
Maximum head elevation	90°
Maximum Trend/Reverse Trend	18°
Floor to base clearance	3.5" (8.9 cm)

Caster size	8" (20.3 cm)
-------------	--------------

Warranty

Frame & Welds for life of product

Three-year parts

Two-year mattress

One-year labor

Accessory Options — Please confirm accessory availability for model of stretcher ordered

Accessory IV pole	Liquid oxygen tank holder
Ankle stirrups	Oxygen tank holder
Armboard	Paper roll dispenser
Chart holder	Patient tray
Convertible footboard	Restraint straps
Detachable footboard	Siderail covers
Foot extender	Stretcher pillow
Infusion support system	Transport straps
IV transporter	Utility tray
Lateral X-ray cassette holder	

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2018 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
202574 rev 2 21-MAY-2018 ENG – US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com

Enhancing outcomes for
patients and their caregivers:

Hill-Rom