





# Hillrom™


HILLROM  
1069 STATE ROUTE 46 EAST  
BATESVILLE, IN 47006  
PHONE: (812) 934-7777  
FAX: (812) 934-8189  
[www.hillrom.com](http://www.hillrom.com)



TOP VIEW



FRONT ELEVATION



SIDE ELEVATION

## SPECIFICATIONS

**ILLUMINATION LEVEL EC AT 1 M (LX):** 80,000

**PATTERN SIZE (D10) AT 1 M (CM) | (INCH):** 16 | 6.3"

**COLOR TEMPERATURE (K):** 4,500


**COLOR RENDERING INDEX RA:** 95

**AVERAGE SERVICE LIFE OF THE LEDS (H):** > 50,000

© [[2018]] HILLROM SERVICES, INC. ALL RIGHTS RESERVED.

## NOTES:

1. INSTALLATION TO BE COMPLETED IN ACCORDANCE WITH MANUFACTURER'S SPECIFICATIONS.
2. DO NOT SCALE DRAWING.
3. THIS DRAWING IS INTENDED FOR USE BY ARCHITECTS, ENGINEERS, CONTRACTORS, CONSULTANTS AND DESIGN PROFESSIONALS FOR PLANNING PURPOSES ONLY. THIS DRAWING MAY NOT BE USED FOR CONSTRUCTION.
4. ALL INFORMATION CONTAINED HEREIN WAS CURRENT AT THE TIME OF DEVELOPMENT BUT MUST BE REVIEWED AND APPROVED BY THE PRODUCT MANUFACTURER TO BE CONSIDERED ACCURATE.
5. CONTRACTOR'S NOTE: FOR PRODUCT AND COMPANY INFORMATION VISIT [www.CADdetails.com/info](http://www.CADdetails.com/info) AND ENTER REFERENCE NUMBER 4794-107



## SURGICAL SOLUTIONS (PREVIOUSLY TRUMPF MEDICAL)

TRULIGHT® 1000 EXAM LIGHT