

Baxter
Welch Allyn
Connex
SPOT MONITOR



MAKE SMARTER, MORE INFORMED CLINICAL DECISIONS.

Meet the **Welch Allyn Connex Spot Monitor**

It's easy to use and supports comprehensive bedside monitoring—from routine vital signs acquisition to early warning scoring and interval monitoring. And we can integrate with your EMR, giving your care teams immediate access to critical information—when and where they need it.



Our guides and features allow you to customize and standardize your workflows.

- Select from six customizable Early Warning Score protocols (eg. MEWS, PEWS, NEWS2) to help detect patient deterioration
- Large color touchscreen and single sign on capability simplifies training and operation
- Manually enter weight, height, respiration rate, pain level and document modifiers, like body position and Is and Os

- Spot-check respiration measurement with the Masimo RRp technology
- Our **SureBP** technology captures BP readings on inflation in about 20-seconds
- Select Masimo or Nellcor
- Choose **Braun ThermoScan® PRO 6000** ear or **SureTemp Plus** oral/axillary/rectal thermometry
- Works with a wide range of cuffs to screen adult, pediatric and neonatal patients

DESIGNED WITH PATIENT OUTCOMES IN MIND.

- 1 Customized features help clinicians identify and respond to signs of patient deterioration sooner
- 2 Vital signs and other patient health data points are sent directly to the EMR to help reduce errors and omissions
- 3 Clinical technologies support and produce readings you can rely on



VISUALLY ASSESSING RESPIRATION RATE? GO DIGITAL.

Digitally spot-check a patient's respiration rate without changing your workflow. Acquire accurate, connected readings at the fingertip helping your team reduce transcription errors and detect signs of deterioration.

SECURE YOUR NETWORK AND PATIENT HEALTH DATA.



If data security is your top concern, we have you covered. When you equip your **Connex** Spot Monitor with **Welch Allyn** EMR connectivity solutions, we align with your network security protocols to help you protect your patients and infrastructure.

- FIPS-approved² encryption of transmitted data with device and server authentication
- Digitally signed software images and configuration files protect devices from unauthorized changes
- Single sign-on compatibility, including Imprivata, supports quick, secure logins
- Screen lock feature protects idle devices from unauthorized access
- Patient data is automatically deleted after it's sent to the EMR and cleared from the device



ACCESSORIES

7000-APM	Accessory Power Management Stand
7000-MS3	Classic Mobile Stand
7000-MWS	Mobile Work Surface Stand
7000-GCX	GCX Wall Channel Mount
6000-NC	VSM 6000 Nurse Call Cable
7000-916HS	2D Barcode Scanner
7000-916HSR	2D Barcode Scanner + High Frequency RFID Reader
BATT22	CSM Lithium-Ion Battery
BATT99	APM Stand Lithium-Ion Battery
4500-35	Blood Pressure Hose with FlexiPort (10 ft)
7000-33	Neonatal Blood Pressure Hose (10 ft)
S1-CSM-PRO-PS	CSM, Service Protection Program, 3-years
S1-CSM-PROPLUS-PS	CSM, Protection+ Program, 3-years

The **Connex** Spot Monitor is compatible with a wide range of cleaning products to help you reduce the risk of cross-contamination.



ACCESSORY POWER MANAGEMENT STAND

7000-APM

Extend your battery life to up to 17 hours of use, helping minimize downtime due to re-charging. Plus, the integrated work surface and built-in accessory/cleaning storage were designed to enhance your productivity.

Want to learn more? Contact your **Baxter** account representative today for a demo or visit hillrom.com/csm.

Rx only. For safe and proper use of the product mentioned herein, please refer to the appropriate Operator’s Manual or Instructions for Use.

References

1. Medical Device and Equipment Risk Assessment. Data on File.
2. Validation OpenSSL FIPS Object Module v2.0. Data on File.

Baxter, Connex, FlexiPort, Hillrom, SureBP, SureTemp Plus, and Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries.

Braun is a trademark of Braun, GmbH.

ThermoScan is a trademark of Helen of Troy Limited.

Any other trademarks, product names or brand images appearing herein are the property of their respective owners.

US-FLC45-230022 (v2.0) 03/2025