




Hillrom™

HILLROM
1069 STATE ROUTE 46 EAST
BATESVILLE, IN 47006
PHONE: (812) 934-7777
FAX: (812) 934-8189
www.hillrom.com

☒ SELECT DESIRED MODEL:

- ☐ 6400 STANDARD OFFERING - 4 HOST USB PORTS FOR ACCESSORIES, NURSE CALL, RADIO READY
- ☐ 6500 WIRELESS - SAME AS 6400 SERIES PLUS ALL MODELS INCLUDE AN INTERNAL 802.11 A/B/G WIRELESS RADIO
- ☐ 6700 STANDARD OFFERING: 4 USB PORTS FOR ACCESSORIES, NURSE CALL, RADIO READY. NIBP AND SPO2 STANDARD. SPOT-CHECK, TRIAGE, INTERVALS MONITORING AND CONTINUOUS MONITORING PROFILE INCLUDED. RESPIRATION OPTION NOT INCLUDED.
- ☐ 6700 STANDARD OFFERING PLUS CO2 - ALL MODELS INCLUDE COVIDIEN(R) CAPNOGRAPHY
- ☐ 6700 STANDARD OFFERING PLUS RRA - ALL MODELS INCLUDE MASIMO ACOUSTIC RESPIRATION
- ☐ 6700 STANDARD OFFERING PLUS EARLYSENSE - ALL MODELS INCLUDE EARLYSENSE CONTACT-FREE MONITORING
- ☐ 6800 STANDARD OFFERING - SAME AS 6700 PLUS ALL MODELS INCLUDE INTERNAL 802.11A/B/G WIRELESS RADIO.
- ☐ 6800 STANDARD OFFERING PLUS CO2 - ALL MODELS INCLUDE COVIDIEN CAPNOGRAPHY
- ☐ 6800 STANDARD OFFERING PLUS RRA - ALL MODELS INCLUDE MASIMO ACOUSTIC RESPIRATION
- ☐ 6800 STANDARD OFFERING PLUS EARLYSENSE - ALL MODELS INCLUDE EARLYSENSE CONTACT-FREE MONITORING




FEATURES:

- SUREBP® SUPPORTS ALL AGES AND PROVIDES 15-SECOND READINGS AS THE CUFF INFLATES.
- SINGLE SIGN-ON OPTIONS AVAILABLE FOR QUICK, SECURE CLINICIAN LOGIN.
- OPTIONAL 3- OR 5-LEAD ECG BEDSIDE MONITORING.
- OPTIONAL COVIDIEN® CAPNOGRAPHY (CO2).
- OPTIONAL MASIMO® ACOUSTIC RESPIRATORY MONITORING (RRA).
- HEALTH O METER®, SECA AND DETECTO WEIGHT SCALE CONNECTIVITY.
- OPTIONAL EARLYSENSE® CONTACT-FREE RESPIRATORY MONITORING AND MOTION SENSING.
- SURETEMP® PLUS AND/OR BRAUN THERMOSCAN® PRO 6000 EAR THERMOMETRY (OPTIONAL).
- MASIMO® OR NELLCOR® PULSE OXIMETRY (SPO2) ALGORITHMS (OPTIONAL).
- OPTIONAL MASIMO® NONINVASIVE CONTINUOUS MONITORING OF TOTAL HEMOGLOBIN (SPHB).
- OPTIONAL CUSTOM DATA INPUTS.
- OPTIONAL EARLY WARNING SCORE INTEGRATION.
- WIRELESS CONNECTIVITY TO LEADING ELECTRONIC MEDICAL RECORDS.
- SHIPS STANDARD WITH ONE ADULT AND ONE LARGE ADULT CUFF.

© [[2018]] HILLROM SERVICES,
INC. ALL RIGHTS RESERVED.

NOTES:

1. INSTALLATION TO BE COMPLETED IN ACCORDANCE WITH MANUFACTURER'S SPECIFICATIONS.
2. DO NOT SCALE DRAWING.
3. THIS DRAWING IS INTENDED FOR USE BY ARCHITECTS, ENGINEERS, CONTRACTORS, CONSULTANTS AND DESIGN PROFESSIONALS FOR PLANNING PURPOSES ONLY. THIS DRAWING MAY NOT BE USED FOR CONSTRUCTION.
4. ALL INFORMATION CONTAINED HEREIN WAS CURRENT AT THE TIME OF DEVELOPMENT BUT MUST BE REVIEWED AND APPROVED BY THE PRODUCT MANUFACTURER TO BE CONSIDERED ACCURATE.
5. CONTRACTOR'S NOTE: FOR PRODUCT AND COMPANY INFORMATION VISIT www.CADdetails.com/info AND ENTER REFERENCE NUMBER 4794-111



PATIENT MONITORING (PREVIOUSLY WELCH ALLYN)

CONNEX® VITAL SIGNS MONITOR