





# Hillrom™


HILLROM  
1069 STATE ROUTE 46 EAST  
BATESVILLE, IN 47006  
PHONE: (812) 934-7777  
FAX: (812) 934-8189  
[www.hillrom.com](http://www.hillrom.com)



TOP VIEW



SIDE ELEVATION



FRONT ELEVATION

## SPECIFICATION

**ILLUMINATION LEVEL EC IN 1M (LX):** 160,000

**PATTERN SIZE (D10) AT 1 M (CM) | (INCH):** 16 – 25 | 6.3" – 9.8"

**PATTERN SIZE CONSTANT OVER WORKING DISTANCE:** STANDARD (31" TO 51")

**COLOR TEMPERATURE (K):** 3,500 / 4,000 / 4,500 / 5,000

**COLOR RENDERING INDEX RA:** MAX. 97

© [[2018]] HILLROM SERVICES, INC. ALL RIGHTS RESERVED.

## NOTES:

1. INSTALLATION TO BE COMPLETED IN ACCORDANCE WITH MANUFACTURER'S SPECIFICATIONS.
2. DO NOT SCALE DRAWING.
3. THIS DRAWING IS INTENDED FOR USE BY ARCHITECTS, ENGINEERS, CONTRACTORS, CONSULTANTS AND DESIGN PROFESSIONALS FOR PLANNING PURPOSES ONLY. THIS DRAWING MAY NOT BE USED FOR CONSTRUCTION.
4. ALL INFORMATION CONTAINED HEREIN WAS CURRENT AT THE TIME OF DEVELOPMENT BUT MUST BE REVIEWED AND APPROVED BY THE PRODUCT MANUFACTURER TO BE CONSIDERED ACCURATE.
5. CONTRACTOR'S NOTE: FOR PRODUCT AND COMPANY INFORMATION VISIT [www.CADdetails.com/info](http://www.CADdetails.com/info) AND ENTER REFERENCE NUMBER 4794-106


**AVERAGE SERVICE LIFE OF THE LEDS (H):** 60,000

**ADAPTIVE LIGHT CONTROL PLUS (ALC PLUS):** STANDARD

**STERILE LIGHT CONTROL (SLC):** OPTIONAL

**ACTIVE LIGHT MANAGEMENT:** STANDARD

**CAMERA PREPARATION:** STANDARD



## SURGICAL SOLUTIONS (PREVIOUSLY TRUMPF MEDICAL)

ILED® 7 SURGICAL LIGHT (SINGLE HEAD)