





Hillrom™

HILLROM
1069 STATE ROUTE 46 EAST
BATESVILLE, IN 47006
PHONE: (812) 934-7777
FAX: (812) 934-8189
www.hillrom.com



FEATURES

- SUPPORTS A WIDE RANGE OF POSITIONS ACROSS SURGICAL APPLICATIONS.
- 1,000LB WEIGHT CAPACITY, COLLISION CONTROL & MORE.
- MODULAR DESIGN, 360° IMAGING & LIGHT MESSAGING SYSTEM.
- SIMPLE REMOTE CONTROL, HOOK COUPLER SYSTEM & MORE.

SPECIFICATIONS

NET BED WEIGHT: 250 KG (551.2 LBS).

HEIGHT ADJUSTMENT: 58.8 CM - 114.8 CM +/- 0.5 CM (23.5" - 45.2" +/- 0.2").

TRENDELENBURG / REVERSE TRENDELENBURG: +/- 35° (+/- 2°).


TILT: +/- 25° (+/- 2°)

LONGITUDINAL SLIDE: 40 CM ± 0.5 CM (15.75" ± 0.2")

HEAD SIDE: 20 CM (7.9")

FOOT SIDE: 20 CM (7.9").


MAXIMUM OVERALL LOAD: 454 KG (1,000 LBS)



© [[2021]] HILLROM SERVICES, INC. ALL RIGHTS RESERVED.

NOTES:

1. INSTALLATION TO BE COMPLETED IN ACCORDANCE WITH MANUFACTURER'S SPECIFICATIONS.
2. DO NOT SCALE DRAWING.
3. THIS DRAWING IS INTENDED FOR USE BY ARCHITECTS, ENGINEERS, CONTRACTORS, CONSULTANTS AND DESIGN PROFESSIONALS FOR PLANNING PURPOSES ONLY. THIS DRAWING MAY NOT BE USED FOR CONSTRUCTION.
4. ALL INFORMATION CONTAINED HEREIN WAS CURRENT AT THE TIME OF DEVELOPMENT BUT MUST BE REVIEWED AND APPROVED BY THE PRODUCT MANUFACTURER TO BE CONSIDERED ACCURATE.
5. CONTRACTOR'S NOTE: FOR PRODUCT AND COMPANY INFORMATION VISIT www.CADdetails.com/info AND ENTER REFERENCE NUMBER 4794-125



SURGICAL SOLUTIONS

PST 500 PRECISION SURGICAL TABLE