




Hillrom™

HILLROM
1069 STATE ROUTE 46 EAST
BATESVILLE, IN 47006
PHONE: (812) 934-7777
FAX: (812) 934-8189
www.hillrom.com

☒ SELECT DESIRED MOUNT:

☐ STANDARD

☐ MOTORIZED BASE AND WIRELESS INTERFACE



SPECIFICATION

MAXIMUM LOAD: 1,000 LBS (450 KG)

HEIGHT ADJUSTMENT: 23 IN – 45 IN (± 0.4 IN) (59 CM – 114 CM (± 1 CM))

LIFT: 22 IN ± 0.4 IN (55 CM ± 1 CM)

TRENDELENBURG / REVERSE TRENDELENBURG: $\pm 45^\circ$

TILT: $\pm 30^\circ$

LONGITUDINAL SLIDE: 18 IN ± 0.4 IN (46 CM ± 1 CM)


ADJUSTMENT RANGE OF LEG SECTION JOINTS: $+ 90^\circ / - 105^\circ$

ADJUSTMENT OF LOWER BACK SECTION JOINTS: $+ 90^\circ / - 45^\circ$

© [[2018]] HILLROM SERVICES, INC. ALL RIGHTS RESERVED.

NOTES:

1. INSTALLATION TO BE COMPLETED IN ACCORDANCE WITH MANUFACTURER'S SPECIFICATIONS.
2. DO NOT SCALE DRAWING.
3. THIS DRAWING IS INTENDED FOR USE BY ARCHITECTS, ENGINEERS, CONTRACTORS, CONSULTANTS AND DESIGN PROFESSIONALS FOR PLANNING PURPOSES ONLY. THIS DRAWING MAY NOT BE USED FOR CONSTRUCTION.
4. ALL INFORMATION CONTAINED HEREIN WAS CURRENT AT THE TIME OF DEVELOPMENT BUT MUST BE REVIEWED AND APPROVED BY THE PRODUCT MANUFACTURER TO BE CONSIDERED ACCURATE.
5. CONTRACTOR'S NOTE: FOR PRODUCT AND COMPANY INFORMATION VISIT www.CADdetails.com/info AND ENTER REFERENCE NUMBER 4794-100



SURGICAL SOLUTIONS (PREVIOUSLY TRUMPF MEDICAL)

TRUSYSTEM® 7000DV SURGICAL TABLE