



CERTIFIED PREVENTIVE MAINTENANCE

Beds and Surfaces

REDUCE SAFETY RISKS AND STAY COMPLIANT WITH:

✓	Documentation <ul style="list-style-type: none"> ▪ Joint commission requirements
✓	Electrical Safety Certification <ul style="list-style-type: none"> ▪ Measure leakage and ground resistance ▪ Inspect power cord and line filter for damage
✓	Functional Controls <ul style="list-style-type: none"> ▪ Ensure caregiver and patient controls are operating correctly
✓	Inspecting Internal Components <ul style="list-style-type: none"> ▪ Air blowers, mattress bladders and hoses ▪ Check mattress interior for moisture or stains
✓	Scale Calibration
✓	Inspecting Safety Features <ul style="list-style-type: none"> ▪ Siderails ▪ CPR feature ▪ Alerts, SideCom® Communication System, Bed Exit, SafeView® + System and wireless connectivity
✓	Inspecting/Adjusting Movement Limits <ul style="list-style-type: none"> ▪ Trendelenburg, Auto Contour™ feature, chair, hi-low, head articulation, foot retraction/articulation, and SlideGuard® feature
✓	Inspecting Bed Movement Functions <ul style="list-style-type: none"> ▪ Brakes, steering and 5th wheel ▪ Powered drive system
✓	Battery Backup <ul style="list-style-type: none"> ▪ Check that battery activates when bed is unplugged and has sufficient power ▪ Confirm caregiver controls function with battery backup
✓	Improving Appearance <ul style="list-style-type: none"> ▪ Paint touch ups ▪ Replace labels and plastic parts

hillrom.com

* The tests listed are typical and customary. All tests may not be applicable to all equipment.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.