


Welch Allyn Connex® Integrated Wall System Mounting Template

REF 105331



Mounting hardware

Mounting diagram

Refer to the *Startup guide* for complete installation instructions.

WARNING Failure to install the mandatory security screw might result in personal injury and equipment damage.

- ① Mounting rail bracket
- ② Wall system
- ③ Braun® thermometer dock (optional)
- ④ Accessory bin bracket
- ⑤ Accessory bin

WelchAllyn®

0 3 6

SCALE: INCHES

Verify scale prior to using template for device installation