


<b>Document Description:</b> Small Platform Stand Release	<b>Document Number:</b> 80025237 <b>Version:</b> A
	
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only.	

## Customer Service Bulletin


<b>Product:</b> <i>Small Platform Universal Stand</i>	<b>Date:</b> 2019-06-13
<b>Subject:</b> Small Platform Stand Release	
<b>HW Version(s) Affected:</b> N/A	<b>SW Version(s) Affected:</b> N/A
<b>Serial Numbers Affected:</b> N/A	<b>Lot or Date Code Affected:</b> N/A

<b>Classification:</b> Informational Only
<b>Distribution:</b> <input checked="" type="checkbox"/> Customer Care <input checked="" type="checkbox"/> Product Service <input checked="" type="checkbox"/> Field Service <input checked="" type="checkbox"/> ASPs <input checked="" type="checkbox"/> Distributors <input checked="" type="checkbox"/> Customers

**Summary:** Hillrom has replaced the SPS2 Stand (Part Number 413020) on configurations with a new small platform stand (Part Number 413171).

**Product Changes:**

- The new small platform stand (Part Number 413171) will have a Welch Allyn “Cool Gray” plastic caster base. The SPS2 Stand (Part Number 413020) had a white cast metal base.
- The new small platform stand (Part Number 413171) will have a plastic pole to pole coupler. The SPS2 Stand (Part Number 413020) had a cast metal pole to pole coupler.
- No Serviceable parts are available for this new stand.
- warranty issues with the previous SPS2 Stand (413020) will use the below service kits or result in a product replacement with the new stand.


<b>Document Description:</b> Small Platform Stand Release	<b>Document Number:</b> 80025237 <b>Version:</b> A
	
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only.	

407431	2" CASTOR SET; 1 BRAKING, 1 NON-BRAKING
708600	MOUNTING PLATE, W/ CAPTIVE SCREW, MS2
406560	MOUNTING OP, SURETEMP PLUS, MOBILE STAND
407645	ASSY, PLATFORM BASKET
407678	PROBP 3400 BASKET ADAPTER WITH SCREW
407679	5W POWER SUPPLY MOUNTING KIT
411125	MOBILE STAND MOUNTING BRACKET, PROBP2400
413020	UNIVERSAL SMALL PLFM MOBILE STAND
411162	MOBILE STAND HARWARE KIT
407867	767 SPS2 MOBILE STAND BRACKET KIT
413063	HOLSTER KIT, PRO-BP 2000


**Mounting/Assembling Changes:**

There are **no** changes in the mounting bracket/device interface.

**PRODUCT IDENTIFIERS:** [Refer to table below](#)

<b>Document Description:</b> Small Platform Stand Release	<b>Document Number:</b> 80025237 <b>Version:</b> A
	
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only.	

LISTED PRODUCT NUMBERS:	Hillrom Part Number	Hillrom Part Description
		2400-MS
	34BFST-2	PROBP 3400, SUREBP, BLUETOOTH, MOBILE, EU
	34BFST-4	PROBP 3400, SUREBP, BLUETOOTH, MOBILE, UK
	34BFST-6	PROBP 3400, SUREBP, BLUETOOTH, MOBILE, AU
	34BFST-7	PROBP 3400, SUREBP, BLUETOOTH, MOBILE, ZA
	34BFST-B	PROBP 3400, SUREBP, BLUETOOTH, MOBILE, NA
	34BFST-C	PROBP 3400, SUREBP, BLUETOOTH, MOBILE, CN
	34BFST-J	PROBP 3400, SUREBP, BLUETOOTH, MOBILE, JP
	34BXST-2	PROBP 3400, BLUETOOTH, MOBILE, EU
	34BXST-4	PROBP 3400, BLUETOOTH, MOBILE, UK
	34BXST-6	PROBP 3400, BLUETOOTH, MOBILE, AU
	34BXST-7	PROBP 3400, BLUETOOTH, MOBILE, ZA
	34BXST-B	PROBP 3400, BLUETOOTH, MOBILE, NA
	34BXST-C	PROBP 3400, BLUETOOTH, MOBILE, CN
	34BXST-J	PROBP 3400, BLUETOOTH, MOBILE, JP
	34XFST-2	PROBP 3400, SUREBP, MOBILE, EU
	34XFST-4	PROBP 3400, SUREBP, MOBILE, UK
	34XFST-6	PROBP 3400, SUREBP, MOBILE, AU
	34XFST-7	PROBP 3400, SUREBP, MOBILE, ZA
	34XFST-B	PROBP 3400, SUREBP, MOBILE, NA
	34XFST-C	PROBP 3400, SUREBP, MOBILE, CN
	34XFST-J	PROBP 3400, SUREBP, MOBILE, JP
	34XXST-2	PROBP 3400, MOBILE, EU
	34XXST-4	PROBP 3400, MOBILE, UK
	34XXST-6	PROBP 3400, MOBILE, AU
	34XXST-7	PROBP 3400, MOBILE, ZA
	34XXST-B	PROBP 3400, MOBILE, NA
	34XXST-C	PROBP 3400, MOBILE, CN
	34XXST-J	PROBP 3400, MOBILE, JP
	4600-61	MOBILE STAND W/MOUNTING KIT PROBP 3400
	7670-03	767 MOBILE GAGE, ADULT, 5 LEG BASE
	7670-03CB	767 MOBILE GAGE, ADULT, 2-PC 5 LEG BASE
	7670-03K	767 MOBILE GAGE, W/O CUFF, 5 LEG BASE
	7670-06P	767 MOBILE STAND KIT - PLATFORM STYLE
	7670-10	767 MOBILE GAGE, ADULT,5 LEG BASE
	7670-10CB	767 MOBILE GAGE, ADULT,2-PC,5 LEG BASE
	7670-10KPA	767 ANEROID (KPA ) ON 5 LEG STAND W/CUFF
	7670-12	MOBILE STAND FOR WA 76710
	107195	STAND KIT, PRO-BP 2000
	01690-700	PKG THERM SYS, M690,9FT OR W/ST
	01692-700	PKG THERM SYS, M692,9FT OR W/ST

<b>Document Description:</b> Small Platform Stand Release	<b>Document Number:</b> 80025237 <b>Version:</b> A
	
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only.	

		08273-000	STAND, M690/M692, W/PROBE&WELL HLDR
		7670-03	767 MOBILE GAGE, ADULT, 5 LEG BASE
		7670-03CB	767 MOBILE GAGE, ADULT, 2-PC 5 LEG BASE
		7670-03K	767 MOBILE GAGE, W/O CUFF, 5 LEG BASE
		7670-10	767 MOBILE GAGE, ADULT,5 LEG BASE
		7670-10CB	767 MOBILE GAGE, ADULT,2-PC,5 LEG BASE
		7670-10KPA	767 ANEROID (KPA ) ON 5 LEG STAND W/CUFF
		7670-10BR	767 MOBILE GAGE, ADULT,5 LEG BASE Brazil

<b>Version</b>	<b>Sec, Pg, Para Changed</b>	<b>Change Made</b>	<b>Date Version Created</b>	<b>Version Created By (initials)</b>
A	N/A	Initial Release	2019-06-13	KMG