



Printed or electronic versions of this document not accessed directly from the designated Hillrom Controlled Quality Information System are For Reference Only.

Customer Service Bulletin

Product:	Connex Spot Monitor (CSM) Connex Vital Signs Monitor (CVSM) Connex Integrated Wall System (CIWS)		
Subject:	Nurse Call Cable 6000-NC Update		
HW Version(s) Affected:	ALL	SW Version(s) Affected:	ALL
Serial Numbers Affected:	ALL	Lot or Date Code Affected:	ALL

Classification:	As Needed		
Distribution:	<input checked="" type="checkbox"/> Customer Care	<input checked="" type="checkbox"/> Product Service	<input checked="" type="checkbox"/> Field Service
	<input checked="" type="checkbox"/> ASPs	<input checked="" type="checkbox"/> Distributors	<input checked="" type="checkbox"/> Customers

Summary:

Due to the obsolescence of the phono plug (part number 406213) used on the end of the VSM 6000 Nurse Call Cable Part number 6000-NC, A new version of this part is being released to support the Vital Signs Monitors (CSM, CVSM, and CIWS). This change has no impact to design functionality. The product IFU's and the instruction card part number 713970 supplied in the box for 6000-NC will be updated.

Current CSM IFU Content:


To connect the monitor to a Nurse Call system, you must have a cable that has been adapted to your Nurse Call system (REF 6000-NC), rated 25V AC or 60V DC maximum at 1A maximum.

Current CVSM/CIWS IFU Content:

To connect the monitor to a Nurse Call system, you must have a cable that has been adapted to your Nurse Call system (REF 6000-NC), rated 50V DC maximum at 500mA (DC or ACRMS).




Current Instruction Card Content for 6000-NC:



CSM, CVSM, and CIWS IFU Updates Specifications for the new cable:

To connect the monitor to a Nurse Call system, you must have a cable that has been adapted to your Nurse Call system (REF 6000-NC), rated 24V AC or 24V DC maximum at 500MA maximum.

Updates for the Instruction Card Content 6000-NC:



END OF BULLETIN

THIS INFORMATION IS THE PROPERTY OF HILLROM HOLDINGS AND AS SUCH SHALL NOT BE REPRODUCED, COPIED, OR USED AS A BASIS FOR THE MANUFACTURE OR SALE OF EQUIPMENT, DEVICES, OR SOFTWARE WITHOUT THE EXPRESS WRITTEN PERMISSION OF HILLROM.



Version	Sec, Pg, Para Changed	Change Made	Date Version Created	Version Created By (initials)
A	N/A	Initial Release	2021-12-07	KMG

THIS INFORMATION IS THE PROPERTY OF HILLROM HOLDINGS AND AS SUCH SHALL NOT BE REPRODUCED, COPIED, OR USED AS A BASIS FOR THE MANUFACTURE OR SALE OF EQUIPMENT, DEVICES, OR SOFTWARE WITHOUT THE EXPRESS WRITTEN PERMISSION OF HILLROM.