

Question or Item	Answer
Did you verify that you are using the most recent version of this form?	☒ Yes ☐ No
Date:	01-AUG-2022
Individual creating the CSB:	Steve Bennett
Product(s) affected:	Diagnostic sets (EU, Canada) 71-PM3LXE & 71-XM3LXE Desk sets (US) 71-PM3LDE-US & 71-XXXLDE
Hardware Version(s) affected:	ALL
Lot or Date Code Affected:	Starting with Lot Code: 22-080 Ending with Lot Code: 22-112

Subject:

Supply chain challenges have affected everyone, even us at Baxter. To help our customers get access to kits with our new PanOptic Plus ophthalmoscope and MacroView Plus otoscope heads in a timely manner, we updated select sets to help insure more timely delivery of these products. The Diagnostic sets for EU & Canada: **71-PM3LXE & 71-XM3LXE** and Desk sets for the US: **71-PM3LDE-US & 71-XXXLDE** normally ship with the 719-3 / LI-ION PLUS USB HANDLE FOR DESK CHARGER. Due to supply-chain issues / long lead times, these kits were *temporarily* updated to ship with the 71900-USB / LI-ION HANDLE WITH USB CHARGER instead of 719-3 / LI-ION PLUS USB HANDLE FOR DESK CHARGER. Kits with the replacement 71900-USB / LI-ION HANDLE WITH USB CHARGER shipped in kits with Lot Codes 22-080 through 22-112.

UPDATE: Diagnostic sets for EU & Canada: 71-PM3LXE & 71-XM3LXE and Desk sets for the US: 71-PM3LDE-US & 71-XXXLDE have returned to shipping with the 719-3 / LI-ION PLUS USB HANDLE FOR DESK CHARGER.

Issue details:

Details:

This temporary change in contents for these materials did not change any end-user SKUs (set part numbers) packaging, labeling, or pricing. Lot codes (detailed above) track which sets have been updated as such.

The 71900-USB is similar to the 719-3 in that they have:

- Same battery
- Same charge times when using the Desk Charger
- Same on-time
- USB charging
- Both are compatible with new 719-DSK and older 71140 desk chargers
- Both are compatible with Welch Allyn 3.5V otoscopes and ophthalmoscopes

The 71900 differs from the 719-3 in that:

- 71900-USB charges more slowly via USB
- 71900-USB has green button on the rheostat that needs to be pressed to twist the rheostat into the ON position
- 71900-USB rheostat turns smoothly from OFF to full brightness. The 719-3 has 5 brightness levels and an OFF position.



Version	Sec, Pg, Para Changed	Change Made	Date Version Created	Version Created By (initials)
A	N/A	Initial Release	2022-03-07	SMB
B	Lot Code Affected / Details	Updated Bulletin with end Lot Code Updated form from previous template to QS20224	2022-08-01	SMB