

**Title:** Welch Allyn Service Tool 1.11.2 Windows 11 Qualification Notes

Question or Item	Answer
Date Created:	2026-05-13
Product(s) affected:	<b>Welch Allyn Service Tool</b>
Hardware Version(s) affected:	N/A
Software Version(s) affected:	1.11.2
Serial Number(s) affected:	N/A
Lot or Date Code Affected:	N/A
<b>Issue details:</b>  Baxter Inc. is pleased to announce the release of the latest Welch Allyn Service Tool (WAST), Version 1.11.2. is fully qualified for Windows 10 and 11 Professional 64-bit editions.	

**Solution/Action:**

This software release addresses the qualification of Welch Allyn Service Tool (WAST) on Windows 10 and 11 operating systems. For Windows 10, drivers are automatically installed by the WAST install package. For Windows 11, customers should download the latest drivers from the NOTE below prior to installing WAST 1.11.2

**PLEASE NOTE:** Install the Latest Drivers **PRIOR** to installing this version of WAST for Windows 11 from <https://www.hillrom.com/en/services/welch-allyn-service-tool/>