

# Welch Allyn Scale-Tronix® Integrated Scales



When integrating a scale into your healthcare facility, you need to trust that it will provide fast, dependable weight measurement that meets your current and future needs.

Scale-Tronix in-floor integrated scales are recessed in the ground, making the platforms level with the floor for easy-access, multi-purpose weighing. The Scale-Tronix 5102 integrated scale displays weight on a stainless steel, wall-mounted digital readout that makes it easy for clinicians to read and accurately record measurements.

Designed to help reduce fall risks for patients while improving clinical workflows—it's what has made **Scale-Tronix a trusted brand for healthcare environments.**



## Durable. Long Lasting. Designed for Demanding Healthcare Environments.



### 5102 Physician's Stand-On Scale

- Low-profile platform with wall-mounted digital readout

### 6102 Flush-Mounted In-Floor Scale

- For wheelchair, chair and stand-on weighing
- Level with floor to help reduce fall risk and improve ease of access
- Bright, easy-to-read display

### 6154 Flush-Mounted In-Floor Scale

- Designed to easily weigh patients on stretchers

## Find the Right Scale for You

	5102 Physician's Stand-On Scale	6102 Flush-Mounted Floor Scale	6154 Flush-Mounted In-Floor Scale
<b>Weighing Capacity</b>	660 lb / 300 kg	825 lb / 375 kg	825 lb / 375 kg
<b>Platform</b>	Anodized aluminum, 18" x 14"	Anodized aluminum, 32" x 36"	Anodized aluminum, 35" x 54"
<b>Optional Accessories</b>	Wall-mount height gauge and power adapter	Printer, thermal printer paper (box, 15 rolls), handrail, battery pack 7.2 V, wall-mount height gauge (stainless steel, IN/CM)	
<b>Readout</b>	Digital LED display in pounds and kilograms, or kilograms only. Wall mounted.	Digital LED display in pounds and kilograms, or kilograms only. Surface or flush mounted.	
<b>Power Source</b>	6 C-size disposable alkaline batteries. Power adapter available.	110 VAC with rechargeable battery backup	
<b>Rough-In Frame Kit</b>	N/A	Allows for a uniform, finished-looking installation	
<b>Printer (Accessory)</b>	N/A	Built-in, paper tape, uses standard thermal paper	
<b>Resolution</b>	0.1 lb / 100 gm		
<b>Automatic Zero</b>	Automatically returns to zero, ready for next patient		
<b>Weight Recall</b>	Recalls last weight		
<b>Reweigh</b>	Recomputes patient's weight while patient is on scale		
<b>Computer Capability</b>	RS-232 output (optional)		

Welch Allyn, Inc.  
4341 State Street Road  
Skaneateles Falls, NY 13153 USA  
(p) 800.535.6663 (f) 315.685.3361