

Welch Allyn Scale-Tronix® Portable Scales



Welch Allyn Scale-Tronix scales are designed for your patients and your demanding healthcare environment.

With low-profile platforms designed to help reduce fall risks, readings that can help clinicians prescribe correct medication dosages with confidence and a wide variety of models, it's easy to see why **Scale-Tronix is a trusted brand for healthcare environments.**



WelchAllyn®

Durable. Long Lasting. Designed for Demanding Healthcare Environments.



5002 Mobile Stand-On Scale

- Versatile and ideal for high-traffic areas
- Low-profile platform designed to help reduce fall risk

5702 Mobile Bariatric Stand-On Scale

- All the benefits of the 5002 model, with a larger platform designed for bariatric patients
- Weighing capacity up to 1,000 lb (454 kg)

6002/6702/6702W/6702SP Wheelchair Scales

- Three scales in one: wheelchair, chair and stand-on
- Easily accommodates most wheelchairs
- Four sizes for increasing capacity

Find the Right Scale for You

	5002 Stand-On Scales	5702 Bariatric Stand-On Scales	6002/6702/6702W/6702SP Wheelchair Scales
Weighing Capacity	880 lb / 400 kg	1,000 lb / 454 kg	880 lb / 400 kg
Platform	Low-profile 18" w x 19.5" d x 1.5" h	Low-profile 24" w x 26" d x 2" h	6002: 24" w x 26" d 6702: 28" w x 30" d 6702W: 28" w x 32" d 6702SP: 35" w x 30" d
Optional Accessories	Handrail (optional in wheelchair models only; standard in others), integrated height gauge, wall-mount height gauge, printer, thermal printer paper, line cord power supply with battery backup		
Resolution	0.1 lb / 100 gm		
Readout	Digital LED display in pounds and kilograms (kilogram-only options available)		
Power Source (Cordless)	6 D-size disposable alkaline batteries. Line cord available.		
Automatic Zero	Automatically returns to zero, ready for next patient		
Weight Recall	Recalls last weight		
Reweigh	Recomputes patient's weight while patient is on scale		
Wheels	Heavy-duty rubber wheels		
Computer Connectivity	RS-232 output (optional)		

Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361