



## 5002 Mobile Stand-On Scale

Durable. Long Lasting. Designed for Demanding Healthcare Environments.

- Versatile and ideal for high-traffic areas
- Low-profile platform and heavy-duty stainless steel handrail designed to help reduce fall risk

### Configurations

<b>5002-XP-B</b>	Standard weight (lb/kg) (X), printer (P) and power adapter with hospital-grade plug (B)
<b>5002-XX-B</b>	Standard weight (lb/kg) (X), data port (X) and power adapter with hospital-grade plug (B)
<b>5002-XX-X</b>	Standard weight (lb/kg) (X), data port (X) and battery power (X)
<b>5002-KP-B</b>	Kg only (K), printer (P) and power adapter with hospital-grade plug (B)
<b>5002-KX-B</b>	Kg only (K), data port (X) and power adapter with hospital-grade plug (B)
<b>5002-KX-X</b>	Kg only (K), data port (X) and battery power (X)

## Specifications

<b>Weighing Capacity</b>	880 lb / 400 kg
<b>Dimensions*</b>	Width: 22" (55.88 cm) Depth: 28" (71.12 cm) Height: 60" (152.4 cm)
<b>Platform</b>	Low profile Width: 18" (45.72 cm) Depth: 19.5" (49.53 cm) Height: 1.5" (3.81 cm)
<b>Resolution</b>	0.1 lb / 100 gm
<b>Readout</b>	Digital LED display in pounds and kilograms (kilogram-only options available)
<b>Power Source (Cordless)</b>	6 D-size disposable alkaline batteries. Line cord available.
<b>Gross Weight*</b>	61 lbs (27.669 kg)
<b>Warranty</b>	One year
<b>Weighing Mechanism</b>	Electronic load cells
<b>Successive Weighing</b>	Allows weighing a continuous stream of patients without zeroing in between
<b>Automatic Zero</b>	Automatically returns to zero, ready for next patient
<b>Weight Recall</b>	Recalls last weight
<b>Reweigh</b>	Recomputes patient's weight while patient is on scale
<b>Computer Connectivity</b>	RS-232 output (optional)
<b>Low Battery Indicator</b>	Indicates battery needs replacing
<b>Audible Signal</b>	Indicates weight reading is entered into memory
<b>Optional Accessories</b>	Integrated height gauge; wall-mount height gauge; printer; thermal printer paper; line cord power supply
<b>Wheels</b>	Heavy-duty rubber wheels
<b>Printer</b>	Built in for models 5002-XP-B and 5002-KP-B. Available as accessory for others. Paper tape, uses standard thermal paper.

\*Weight and dimensions represent approximate values for box weight and dimensions.

Contact us at [scale-tronix@welchallyn.com](mailto:scale-tronix@welchallyn.com)  
or visit [www.welchallyn.com](http://www.welchallyn.com) to learn more.

Welch Allyn, Inc.  
4341 State Street Road  
Skaneateles Falls, NY 13153 USA  
(p) 800.535.6663 (f) 315.685.3361

[WWW.WELCHALLYN.COM](http://WWW.WELCHALLYN.COM)

© 2017 Welch Allyn MC13742 80021052 Ver A 2017-04

**WelchAllyn®**