



5702 Mobile Bariatric Stand-On Scale

Durable. Long Lasting. Designed for Demanding Healthcare Environments.

- All the benefits of the 5002 model, with a larger platform designed for bariatric patients
- Weighing capacity up to 1,000 lbs (454 kg)
- Versatile and ideal for high-traffic areas
- Low-profile platform and heavy-duty stainless steel handrail designed to help reduce fall risk

Configurations

5702-KP-B	Kg only (K), printer (P) and power adapter with hospital-grade plug (B)
5702-KX-B	Kg only (K), data port (X) and power adapter with hospital-grade plug (B)
5702-KX-X	Kg only (K), data port (X) and battery power (X)
5702-XP-B	Standard weight (lb/kg) (X), printer (P) and power adapter with hospital-grade plug (B)
5702-XX-B	Standard weight (lb/kg) (X), data port (X) and power adapter with hospital-grade plug (B)
5702-XX-X	Standard weight (lb/kg) (X), data port (X) and battery power (X)

Specifications

Weighing Capacity	1,000 lb / 454 kg
Dimensions*	Width: 25" (63.5 cm) Depth: 20" (50.8 cm) Height: 53" (134.62 cm)
Platform	Low-profile Width: 24" (60.96 cm) Depth: 26" (66.04 cm) Height: 2" (5.08 cm)
Resolution	0.1 lb / 100 gm
Readout	Digital LED display in pounds and kilograms (kilogram-only options available)
Power Source (Cordless)	6 D-size disposable alkaline batteries. Line cord available.
Gross Weight*	92 lbs (41.731 kg)
Warranty	One year
Weighing Mechanism	Electronic load cells
Successive Weighing	Allows weighing a continuous stream of patients without zeroing in between
Automatic Zero	Automatically returns to zero, ready for next patient
Weight Recall	Recalls last weight
Reweigh	Recomputes patient's weight while patient is on scale
Computer Connectivity	RS-232 output (optional)
Audible Signal	Indicates weight reading is entered into memory
Low Battery Indicator	Indicates battery needs replacing
Optional Accessories	Integrated height gauge (cm and inches); wall-mount height gauge; printer; thermal printer paper; line cord power supply
Wheels	Heavy-duty rubber wheels
Printer	Built in for models 5702-KP-B and 5702-XP-B. Available as accessory for others. Paper tape, uses standard thermal paper.

*Weight and dimensions represent approximate values for box weight and dimensions.

Contact us at scale-tronix@welchallyn.com
or visit www.welchallyn.com to learn more.

Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361

WWW.WELCHALLYN.COM

© 2017 Welch Allyn MC13743 80021053 Ver A 2017-04

WelchAllyn®