



## 6154 Flush-Mounted In-Floor Scale

Durable. Long Lasting. Designed for Demanding Healthcare Environments.

- Designed to easily weigh patients on stretchers

### Configurations

<b>6154-XPE-B</b>	Standard weight (lb/kg) (X), printer (P) and external readout (E)
<b>6154-XXE-B</b>	Standard weight (lb/kg) (X), data port (X) and external readout (E)
<b>6154-XPF-B</b>	Standard weight (lb/kg) (X), printer (P) and flush-mounted readout (F)
<b>6154-XXF-B</b>	Standard weight (lb/kg) (X), data port (X) and flush-mounted readout (F)
<b>6154-KPE-B</b>	Kg only (K), printer (P) and external readout (E)
<b>6154-KXE-B</b>	Kg only (K), data port (X) and external readout (E)
<b>6154-KPF-B</b>	Kg only (K), printer (P) and flush-mounted readout (F)
<b>6154-KXF-B</b>	Kg only (K), data port (X) and flush-mounted readout (F)

## Specifications

<b>Weighing Capacity</b>	825 lb / 375 kg
<b>Dimensions*</b>	Width: 37" (93.98 cm) Height: 56" (142.24 cm) Depth: 20" (50.8 cm)
<b>Platform</b>	Anodized aluminum 35" (88.9 cm) x 54" (137.16 cm)
<b>Resolution</b>	0.1 lb / 100 gm
<b>Readout</b>	Digital LED display in pounds and kilograms, or kilograms only. Surface or flush mounted.
<b>Power Source</b>	110 VAC with rechargeable battery backup
<b>Gross Weight*</b>	124 lb (56.246 kg)
<b>Rough-In Frame Kit</b>	Allows for a uniform, finished-looking installation
<b>Warranty</b>	One year
<b>Weighing Mechanism</b>	Electronic load cells
<b>Successive Weighing</b>	Allows weighing a continuous stream of patients without zeroing in between
<b>Automatic Zero</b>	Automatically returns to zero, ready for next patient
<b>Weight Recall</b>	Recalls last weight
<b>Reweigh</b>	Recomputes patient's weight while patient is on scale
<b>Computer Connectivity</b>	RS-232 output (optional)
<b>Audible Signal</b>	Indicates weight reading is entered into memory
<b>Optional Accessories</b>	Printer, thermal printer paper (box, 15 rolls), handrail, battery pack 7.2 V, wall-mount height gauge (stainless steel, in/cm)
<b>Printer</b>	Built in for models 6154-XPE-B, 6154-XPF-B, 6154-KPE-B and 6154-KPF-B. Available as accessory for others. Paper tape, uses standard thermal paper.

\*Weight and dimensions represent approximate values for box weight and dimensions.

Contact us at [scale-tronix@welchallyn.com](mailto:scale-tronix@welchallyn.com)  
or visit [www.welchallyn.com](http://www.welchallyn.com) to learn more.

Welch Allyn, Inc.  
4341 State Street Road  
Skaneateles Falls, NY 13153 USA  
(p) 800.535.6663 (f) 315.685.3361

[WWW.WELCHALLYN.COM](http://WWW.WELCHALLYN.COM)

© 2017 Welch Allyn MC13745 80021055 Ver A 2017-04

**WelchAllyn**<sup>®</sup>