



5202 Stow-A-Weigh® Stand-On Scale

Durable. Long Lasting. Designed for Demanding Healthcare Environments.

- Mounts snug to wall with slender 4” profile when not in use
- Weighing platform is counter-balanced to fold down easily—no bending or lifting required
- Low-profile platform helps reduce patient fall risk

Configurations

5202-X	Standard weight (lb/kg) (X)
5202-K	Kg only (K)

Specifications

Weighing Capacity	660 lb / 300 kg
Dimensions*	<p>Platform Width: 18" (45.72 cm) Depth: 18" (45.72 cm) Height: 1.5" (3.81 cm)</p> <p>Overall Width: 27" (68.58 cm) Depth: 4" (10.16 cm) Height: 27" (68.58 cm)</p>
Platform	Low profile—folds against the wall. Counterbalanced for fingertip operation.
Resolution	0.1 lb / 100 gm
Readout	Digital LED display in pounds and kilograms, or kilograms only
Power Source (Cordless)	6 D-sized disposable alkaline batteries
Gross Weight*	46 lb (20.865 kg)
Warranty	One year
Weighing Mechanism	Electronic load cells
Successive Weighing	Allows weighing a continuous stream of patients without zeroing in between
Automatic Zero	Automatically returns to zero, ready for next patient
Weight Recall	Recalls last weight
Reweigh	Recomputes patient's weight while patient is on scale
Computer Connectivity	N/A
Low Battery Indicator	Indicates battery needs replacing
Optional Accessories	Height gauge (wall mount, stainless steel, in/cm)
Height Gauge (Accessory)	Wall mount only

*Weight and dimensions represent approximate values for box weight and dimensions.

Contact us at scale-tronix@welchallyn.com
or visit www.welchallyn.com to learn more.

Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361

WWW.WELCHALLYN.COM

© 2017 Welch Allyn MC13748 80021056 Ver A 2017-04

WelchAllyn®