

Connex® Spot Monitor Mounting and Mobility Options

Vital signs data within your reach—across care settings and workflows

Mobile stands support clinical workflows that require a shared device across multiple rooms.



	Standard Basket Stand Simple portability with a small footprint	Work Surface and Accessory Management Stand Organized storage and cable management, plus a compact mobile work surface footprint	Accessory Power Management Stand Mobile work surface and integrated battery to power through a busy day
--	--	---	--

Part Number*	7000-MS3	7000-MWS	7000-APM
Height	49" (124.46 cm) total 38" (96.52 cm) to top of basket handle	47" (119.38 cm) total 39.25" (99.7 cm) to work surface	47" (119.38 cm) total 41" (104.14 cm) to work surface
Base Diameter	23" (58.42 cm)	23" (58.42 cm)	23" (58.42 cm)
Work Surface Area	N/A	11.62" (29.51 cm) wide 11.25" (28.56 cm) deep	11.5" (29.21 cm) wide 10" (25.4 cm) deep
Accessory Basket	12.5" (31.75 cm) wide 8.25" (20.96 cm) deep 7.25" (18.42 cm) high	11.75" (29.85 cm) wide 8.25" (20.96 cm) deep 7.5" high (19.05 cm) 4 accessory bays	11.75" (29.85 cm) wide 8.25" (20.96 cm) deep 7.5" high (19.05 cm) 4 accessory bays
Weight	13.8 kg / 30.5 lb (fully loaded with device and accessories)	15.4 kg / 34 lb (fully loaded with device and accessories)	18.4 kg / 40 lb (fully loaded with device and accessories)
Power Accessory Functionality	N/A	N/A	9 cell Lithium Ion Battery (Field replaceable, integrated power supply, integrated work light, battery life indicator display, two USB ports)

*Connex Spot Monitor sold separately.

Wall mount options are ideal for care settings that need a dedicated device in every exam room, while desktop stands can be moved as needed.



Desktop Stand

Compact and portable tabletop mount suitable for placement on any work surface

Wall Mount with Basket

Easy-to-access storage and cable management

Integrated Wall System

Complete set of physical exam tools in one integrated package

Part Number(s)*	7000-DST	7000-GCX	77794-MCSM**, 77794-2MCSM***
Dimensions	Height 9.5" (24.13 cm) Width 7.25" (18.42 cm) Depth 4.75" (12.07 cm)	Wall Channel: 19" (48.26 cm) Pivot Arm: 8" (20.32 cm) Wall Mount Utility Basket: 12.75" (32.39 cm) wide x 7.5" (19.05 cm) deep x 5.9" (14.99 cm) high	Height: 12" (30.48 cm) Length: 34" (86.36) Depth*: 2" (5.08 cm) *Depth dimension is for wall board only

The Connex Spot Monitor uses a standard VESA mount. Aftermarket mounting options are available.



Learn more at www.welchallyn.com.

* Connex Spot Monitor sold separately.

** 77794-MCSM model includes mounting bracket for Connex Spot Monitor, Coaxial Ophthalmoscope (#11720), MacroView™ Otoscope (#23810), KleenSpec® Disposable Specula dispenser, vital signs accessory basket, vital signs monitor power supply (#7000-PS) and wall mount panel.

*** 77794-2MCSM model includes mounting bracket for Connex Spot Monitor, PanOptic™ Ophthalmoscope with cobalt-blue filter and corneal viewing lens (#11820), MacroView Otoscope (#23810), KleenSpec Disposable Specula dispenser, vital signs accessory basket, vital signs monitor power supply (#7000-PS) and wall mount panel.