



Hillrom™

Welch Allyn® Multi-Parameter Vital Signs Devices



Spot Vital Signs® 4400 Device



Connex® Spot® Monitor



Connex Vital Signs Monitor



Connex Integrated Wall System

| | | | | |
|---------------------------------------|---|--|--|---|
| Workflows | Spot | Spot, Office, Intervals | Spot, Office, Intervals, Continuous | Spot, Office, Intervals, Continuous |
| Device Type | Basic | Intermediate | Advanced | Advanced |
| Patient Populations | Adult, Pediatric | Adult, Pediatric, Neonatal | Adult, Pediatric, Neonatal | Adult, Pediatric, Neonatal |
| Device Alarms | — | Yes | Yes | Yes |
| Third-Party Alarm Notification | — | — | Yes | Yes |
| Parameters | Blood Pressure, Temperature, SpO ₂ | Blood Pressure, Temperature, SpO ₂ | Blood Pressure, Temperature, SpO ₂ , etCO ₂ , SpHb, ECG, RRa | Blood Pressure, Temperature, SpO ₂ , SpHb, ECG, Otoscope, Ophthalmoscope |
| SureBP® Technology (15-sec) | Yes | Yes | Yes | Yes |
| Blood Pressure Averaging | Yes | Yes | Yes | Yes |
| Diagnostic Handles | — | — | — | Otoscope, Ophthalmoscope |
| SpO₂ Options | Nonin® | Nellcor®, Masimo®, Nonin | Nellcor, Masimo | Nellcor, Masimo |
| Temperature Options | SureTemp® Technology | SureTemp Plus, Braun PRO 6000 | SureTemp Plus, Braun PRO 6000 | SureTemp Plus, Braun PRO 6000 |
| Early Warning Scoring | — | Up to 6 protocols | Up to 6 protocols | Up to 6 protocols |
| Custom Documentation Fields | 5 Fixed and Selectable | Up to 20 | Up to 20 | Up to 20 |
| Wireless Connectivity | — | WiFi®, Bluetooth® | WiFi | WiFi |
| Wired Connectivity | USB | USB, Ethernet | USB, Ethernet | USB, Ethernet |
| Display | 7" Color LCD, Touchscreen | 7" Color LCD, Touchscreen | 8.9" Color LCD, Touchscreen | 8.9" Color LCD, Touchscreen |
| Scale Integration | — | Yes | Yes | Yes |
| Mounting Options | Desktop Stand, Mobile Stand | Power Management Stand, Classic Stand, Work Surface Stand, Wall Mount, Integrated Wall System, Desktop Stand | Cable Management Stand, Basic Stand, Wall Channel | Wall Mounted Unit |
| Accessory Options | — | Scales, Barcode, RFID | Scales, Barcode | Scales, Barcode |
| Memory (# Readings) | 50 | 400 | 400 | 400 |

For more information, please contact your Hillrom representative at 1-800-535-6663.

hillrom.com

Hillrom reserves the right to make changes without notice in design, specifications and models. The only warranty Hillrom makes is the express written warranty extended on the sale or rental of its products.