



WAM™ 12-lead ECG Wireless Acquisition Module

FEATURE	DESCRIPTION
Transmission Protocol	Secure, bidirectional frequency-hopping, beacon and response method links each WAM module to a single electrocardiograph or stress system
Sampling Rate	40,000 samples/second/channel acquisition, 1,000 samples/second/channel transmitted for analysis
User Interface	Three-button operation: On/Off, 12-lead ECG acquisition and rhythm strip acquisition
Battery	1 AA alkaline
Size	4.45 x 4.25 x 1.1" (11.3 x 10.8 x 2.79 cm.)
Weight	6.7 oz. (190 g) with battery

Learn more at www.welchallyn.com.

WelchAllyn®

© 2018 Welch Allyn MC14717 80022563 Ver A 2018-01