



Welch Allyn
Scale-Tronix® 4802D Pediatric/Infant Scale



In a demanding health care environment, depend on the **Scale-Tronix 4802D Pediatric/Infant Scale**, designed and built tough to exceed your expectations.

- 1 Durable**
Heavy-duty stainless-steel construction
- 2 Comprehensive**
Supports daily and serial weighing of neonates and infants
- 3 Reliable**
Memory recall will display the last weight taken
- 4 Versatile**
Automatically converts to one gram resolution below one kilogram to facilitate diaper weighing



User-friendly

The default resolution is user-adjustable.

Easily switch between imperial and metric units or choose a configuration which only displays metric units.



Pounds



Pounds and Ounces



Kilograms



Ounces

Pick the cradle that's most appropriate for your needs.

Standard Cradle



Oversized Cradle



Four-sided Cradle



Order additional cradles separately; they're easy for anyone on staff to swap out as your needs change.

Optional mobile cart allows you to easily relocate your scale

Side walls secure and stabilize scale onto cart platform

Convenient storage for scale liners

Large swivel castors and locking wheels support mobility and stabilization



Scale-Tronix 4802D Pediatric/Infant Scale Product Specifications

Weighing Range	10 gm to 20 kg (0.35 oz to 44 lbs)		
Product Dimensions	Interchangeable Acrylic Cradles: Standard: 25" w x 14" d Oversized: 32" w x 14" d Four-sided: 25" w x 14" d	Stainless Steel Scale Body: Height: 7" Width: 27.5" Depth: 16"	Optional Mobile Cart: Cart Height: 39.5" h Cart Base : 23.5" w x 23.5" d Cart Shelf: 33.5" w x 17.5" d
Display Resolution	<1kg: .04oz/1 gm ≥ 1kg: 0.2 oz/5 gm (default) User Adjustable 1, 2, 5, or 10 gm		
Display Units	Imperial and Metric: lb, lb-oz, oz, kg Metric Only: kg		
Display Readout Type	Green 7-Segment LED		
Power Source	Battery and AC powered 6 D-sized alkaline batteries and power supply included		
Automatic Zero	Automatically returns to zero, ready for next patient		
Recall	Recalls last weight		
Reweigh	Remeasures patient's weight while patient is on scale		
Data Transfer	USB Port (Type B) - Welch Allyn Communication Protocol to Welch Allyn Connex® Vital Signs Device (CSM/CVSM) - ASCII Serial Protocol to Windows 10 environment		

Ordering Information:

STANDARD CRADLE

4802D-AX-XB 4802D Pediatric Scale with Standard Cradle, Imperial and Metric (lb/oz/kg)

4802D-AK-XB 4802D Pediatric Scale with Standard Cradle, Metric (kg) only

OVERSIZED CRADLE

4802D-AX-VB 4802D Pediatric Scale with Oversized Cradle, Imperial and Metric (lb/oz/kg)

4802D-AK-VB 4802D Pediatric Scale with Oversized Cradle, Metric (kg) only

FOUR-SIDED CRADLE

4802D-AX-FB 4802D Pediatric Scale with Four-sided Cradle, Imperial and Metric (lb/oz/kg)

4802D-AK-FB 4802D Pediatric Scale with Four-sided Cradle, Metric (kg) only

ADDITIONAL ITEMS AND ACCESSORIES

48217 Scale Liners (500/case)

412490 Mobile Cart

726551 Power Supply

48200W Standard Cradle only

48222 Oversized Cradle only

48220W Four-sided Cradle only

Contact us at scale-tronix@welchallyn.com or visit www.welchallyn.com to learn more.



Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361